



# GPR62 (C-16): sc-54816

## BACKGROUND

GPR62 (probable G protein-coupled receptor 62, hGPCR8) is a 368 amino acid protein encoded by the human GPR62 gene. GPR62 is an orphan receptor member of the G protein-coupled receptor 1 family. G protein-coupled receptors (GPCRs or GPRs) contain seven transmembrane domains and transduce extracellular signals through heterotrimeric G proteins. Key roles for G protein-coupled receptors include control of protein maturation and cell surface delivery and providing the correct framework for interactions with both heterotrimeric G proteins and arrestins to allow signal generation and its termination. GPR62 is expressed in brain tissue, most notably the basal forebrain, frontal cortex, caudate, putamen, thalamus and hippocampus.

## REFERENCES

- Lee, D.K., George, S.R., Cheng, R., Nguyen, T., Liu, Y., Brown, M., Lynch, K.R. and O'Dowd, B.F. 2001. Identification of four novel human G protein-coupled receptors expressed in the brain. *Brain Res. Mol. Brain Res.* 86: 13-22.
- Cikos, S., Gregor, P. and Koppel, J. 2001. Cloning of a novel biogenic amine receptor-like G protein-coupled receptor expressed in human brain. *Biochim. Biophys. Acta* 1521: 66-72.
- Takeda, S., Kadowaki, S., Haga, T., Takaesu, H. and Mitaku, S. 2002. Identification of G protein-coupled receptor genes from the human genome sequence. *FEBS Lett.* 520: 97-101.
- Conner, A.C., Hay, D.L., Simms, J., Howitt, S.G., Schindler, M., Smith, D.M., Wheatley, M. and Poyner, D.R. 2004. A key role for transmembrane prolines in calcitonin receptor-like receptor agonist binding and signalling: implications for family B G protein-coupled receptors. *Mol. Pharmacol.* 67: 20-31.
- Gregory, S.G., Barlow, K.F., McLay, K.E., Kaul, R., Swarbreck, D., Dunham, A., Scott, C.E., Howe, K.L., Woodfine, K., Spencer, C.C., Jones, M.C., Gillson, C., Burrill, W.D., Clegg, S.M., Dhami, P., Dovey, O., Faulkner, L.M., Gribble, S.M., Langford, C.F., Pandian, R.D., Porter, K.M. and Prigmore, E. 2006. The DNA sequence and biological annotation of human chromosome 1. *Nature* 441: 315-321.
- Milligan, G. 2007. A day in the life of a G protein-coupled receptor: the contribution to function of G protein-coupled receptor dimerization. *Br. J. Pharmacol.* 153 (Suppl. 1): S216-229.
- Oldham, W.M. and Hamm, H.E. 2007. Heterotrimeric G protein activation by G protein-coupled receptors. *Nat. Rev. Mol. Cell Biol.* 9: 60-71.

## CHROMOSOMAL LOCATION

Genetic locus: GPR62 (human) mapping to 3p21.1.

## SOURCE

GPR62 (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a C-terminal cytoplasmic domain of GPR62 of human origin.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-54816 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

GPR62 (C-16) is recommended for detection of GPR62 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for GPR62 siRNA (h): sc-62399.

Molecular Weight of GPR62: 37 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.