



MOAP1 (D-14): sc-54850

BACKGROUND

MOAP1 (Modulator of apoptosis 1) is a 352 amino acid protein encoded by the human gene MOAP1. MOAP1 belongs to the PNMA family and contains one BH3-like domain and one RASSF1-binding domain. It is required for death receptor-dependent apoptosis. When MOAP1 is associated with RASSF1, it promotes a Bax conformational change and translocation to mitochondrial membranes in response to TNF and TNFSF10 stimulation. MOAP1 is a homodimer and under normal circumstances, held in an inactive conformation by an intramolecular interaction. Binding to RASSF1 isoform A (RASSF1A) relieves this inhibitory interaction and allows further binding to Bax. It will also bind to BCL2 and BCLX.

REFERENCES

1. Tan, K.O., Tan, K.M., Chan, S.L., Yee, K.S., Bevort, M., Ang, K.C. and Yu, V.C. 2001. MAP-1, a novel proapoptotic protein containing a BH3-like motif that associates with Bax through its Bcl-2 homology domains. *J. Biol. Chem.* 276: 2802-2807.
2. Tan, K.O., Fu, N.Y., Sukumaran, S.K., Chan, S.L., Kang, J.H., Poon, K.L., Chen, B.S. and Yu, V.C. 2005. MAP-1 is a mitochondrial effector of Bax. *Proc. Natl. Acad. Sci. USA* 102: 14623-14628.
3. Baksh, S., Tommasi, S., Fenton, S., Yu, V.C., Martins, L.M., Pfeifer, G.P., Latif, F., Downward, J. and Neel, B.G. 2005. The tumor suppressor RASSF1A and MAP-1 link death receptor signaling to Bax conformational change and cell death. *Mol. Cell* 18: 637-650.
4. Tretyakova, I., Zolotukhin, A.S., Tan, W., Bear, J., Propst, F., Ruthel, G. and Felber, B.K. 2005. Nuclear export factor family protein participates in cytoplasmic mRNA trafficking. *J. Biol. Chem.* 280: 31981-31990.
5. Schüller, M., Jenne, D. and Voltz, R. 2005. The human PNMA family: novel neuronal proteins implicated in paraneoplastic neurological disease. *J. Neuroimmunol.* 169: 172-176.
6. Vos, M.D., Dallol, A., Eckfeld, K., Allen, N.P., Donninger, H., Hesson, L.B., Calvisi, D., Latif, F. and Clark, G.J. 2006. The RASSF1A tumor suppressor activates Bax via MOAP1. *J. Biol. Chem.* 281: 4557-4563.
7. Fu, N.Y., Sukumaran, S.K. and Yu, V.C. 2007. Inhibition of ubiquitin-mediated degradation of MOAP1 by apoptotic stimuli promotes Bax function in mitochondria. *Proc. Natl. Acad. Sci. USA* 104: 10051-10056.

CHROMOSOMAL LOCATION

Genetic locus: MOAP1 (human) mapping to 14q32; Moap1 (mouse) mapping to 12 E.

SOURCE

MOAP1 (D-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MOAP1 of mouse origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-54850 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MOAP1 (D-14) is recommended for detection of MOAP1 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MOAP1 siRNA (m): sc-62630.

Molecular Weight of MOAP1: 40 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.