

Aminoacylase-1 (T-15): sc-55061

BACKGROUND

Aminoacylase-1, also designated N-acyl-L-amino-acid amidohydrolase or ACY-1, is a member of the largest metallopeptidase family called M20A. Aminoacylase-1 is a zinc-binding homodimeric enzyme expressed in kidney, brain, placenta and spleen. It is the most abundant of the aminoacylases. Aminoacylase-1 localizes to the cytoplasm and catalyzes the hydrolysis of N-acylated or N-acetylated amino acids. In addition, Aminoacylase-1 is a sphingosine kinase 1 (SphK1)-interacting protein and may also play a role in regulating responses of the cell to oxidative stress. The gene encoding Aminoacylase-1 is evolutionarily conserved in fish, frog, mouse, rat and human. Deficiency in the Aminoacylase-1 protein may result in defects of brain metabolism and function.

REFERENCES

- Cook, R.M., Burke, B.J., Buchhagen, D.L., Minna, J.D. and Miller, Y.E. 1993. Human Aminoacylase-1. Cloning, sequence, and expression analysis of a chromosome 3p21 gene inactivated in small cell lung cancer. *J. Biol. Chem.* 268: 17010-17017.
- Lindner, H., Berens, W., Kraus, I. and Röhm, K.H. 2001. Mutational analysis of two PWW sequence motifs in human aminoacylase 1. *Biol. Chem.* 381: 1055-1061.
- Lindner, H.A., Lunin, V.V., Alary, A., Hecker, R., Cygler, M. and Menard, R. 2003. Essential roles of zinc ligation and enzyme dimerization for catalysis in the aminoacylase-1/M20 family. *J. Biol. Chem.* 278: 44496-44504.
- Maceyka, M., Nava, V.E., Milstien, S. and Spiegel, S. 2004. Aminoacylase-1 is a sphingosine kinase 1-interacting protein. *FEBS Lett.* 568: 30-34.
- Lindner, H.A., Alary, A., Boju, L.I., Sulea, T. and Menard, R. 2005. Roles of dimerization domain residues in binding and catalysis by Aminoacylase-1. *Biochemistry* 44: 15645-15651.
- Liu, Z., Zhen, Z., Zuo, Z., Wu, Y., Liu, A., Yi, Q. and Li, W. 2006. Probing the catalytic center of porcine Aminoacylase-1 by site-directed mutagenesis, homology modeling and substrate docking. *J. Biochem.* 139: 421-430.
- Sass, J.O., Mohr, V., Olbrich, H., Engelke, U., Horvath, J., Fliegau, M., Loges, N.T., Schweitzer-Krantz, S., Moebus, R., Weiler, P., Kispert, A., Superti-Furga, A., Wevers, R.A. and Omran, H. 2006. Mutations in ACY, the gene encoding Aminoacylase-1, cause a novel inborn error of metabolism. *Am. J. Hum. Genet.* 78: 401-409.
- Engelke, U.F., Sass, J.O., Van Coster, R.N., Gerlo, E., Olbrich, H., Krywawych, S., Calvin, J., Hart, C., Omran, H. and Wevers, R.A. 2007. NMR spectroscopy of Aminoacylase-1 deficiency, a novel inborn error of metabolism. *NMR Biomed.* 21: 138-147.

CHROMOSOMAL LOCATION

Genetic locus: *Acy1* (rat) mapping to 8q32.

SOURCE

Aminoacylase-1 (T-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Aminoacylase-1 of rat origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-55061 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Aminoacylase-1 (T-15) is recommended for detection of Aminoacylase-1 of rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight (predicted) of Aminoacylase-1: 45 kDa.


Molecular Weight (observed) of Aminoacylase-1: 42 kDa.

Positive Controls: rat brain extract: sc-2392.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Aminoacylase-1 (T-15): sc-55061. Western blot analysis of Aminoacylase-1 expression in rat brain tissue extract.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.