

# MPHOSPH6 (S-19): sc-55134

## BACKGROUND

Progression of cells from interphase to mitosis involves alterations in cell structures and activities. The transition from G<sub>2</sub> to M phase is induced by M phase-promoting factor (MPF). In M phase, many proteins are phosphorylated directly by MPF or indirectly by kinases activated by MPF. These M phase phosphoproteins (MPPs), also known as MPHOSPHs, permit disassembly of interphase structures and generation of M phase enzymatic activities and structures. MPP6 (M phase phosphoprotein 6) is a 160 amino acid protein encoded by the human gene MPHOSP6. MPP6, a member of the MPP family, contains one nuclear localization signal motif.

## REFERENCES

1. Matsumoto-Taniura, N., Pirollet, F., Monroe, R., Gerace, L. and Westendorf, J.M. 1997. Identification of novel M phase phosphoproteins by expression cloning. *Mol. Biol. Cell* 7: 1455-1469.
2. Chen, C.Y., Gherzi, R., Ong, S.E., Chan, E.L., Raijmakers, R., Pruijn, G.J., Stoecklin, G., Moroni, C., Mann, M. and Karin, M. 2001. AU binding proteins recruit the exosome to degrade ARE-containing mRNAs. *Cell* 107: 451-464.
3. Leonoudakis, D., Conti, L.R., Anderson, S., Radeke, C.M., McGuire, L.M., Adams, M.E., Froehner, S.C. and Vandenberg, C.A. 2004. Protein trafficking and anchoring complexes revealed by proteomic analysis of inward rectifier potassium channel (Kir2.x)-associated proteins. *J. Biol. Chem.* 279: 22331-22346.
4. Lehner, B. and Sanderson, C.M. 2004. A protein interaction framework for human mRNA degradation. *Genome Res.* 14: 1315-1323.
5. Schilders, G., Raijmakers, R., Raats, J.M. and Pruijn, G.J. 2006. MPP6 is an exosome-associated RNA-binding protein involved in 5.8S rRNA maturation. *Nucleic Acids Res.* 33: 6795-6804.
6. Schilders, G., van Dijk, E. and Pruijn, G.J. 2007. C1D and hMtr4p associate with the human exosome subunit PM/ScI-100 and are involved in pre-rRNA processing. *Nucleic Acids Res.* 35: 2564-2572.

## CHROMOSOMAL LOCATION

Genetic locus: MPHOSPH6 (human) mapping to 16q23.3; Mphosph6 (mouse) mapping to 8 E1.

## SOURCE

MPHOSPH6 (S-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of MPHOSPH6 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-55134 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## APPLICATIONS

MPHOSPH6 (S-19) is recommended for detection of MPHOSPH6 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MPHOSPH6 (S-19) is also recommended for detection of MPHOSPH6 in additional species, including equine, bovine and porcine.

Suitable for use as control antibody for MPHOSPH6 siRNA (h): sc-62635, MPHOSPH6 siRNA (m): sc-62636, MPHOSPH6 shRNA Plasmid (h): sc-62635-SH, MPHOSPH6 shRNA Plasmid (m): sc-62636-SH, MPHOSPH6 shRNA (h) Lentiviral Particles: sc-62635-V and MPHOSPH6 shRNA (m) Lentiviral Particles: sc-62636-V.

Molecular Weight of MPHOSPH6: 19 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.