

TSC-22 D2 (C-20): sc-55406

BACKGROUND

Transforming growth factor β -stimulated clone-22 (TSC-22) acts as a transcriptional regulator to modulate cell growth and differentiation as well as cell death. TSC-22 contains a leucine zipper domain as well as a nuclear export signal, resulting in cytoplasmic localization in living cells. However, concomitant with the induction of apoptosis, TSC-22 translocates from the cytoplasm to the nucleus and shows transcriptional regulatory activity. TSC-22 acts as a major downstream component in both the TGF β pathway and the PPAR γ signaling pathway. The association of these two pathways with tumor suppression and the significant downregulation of TSC-22 mRNA in various cancer types implies an antiproliferative role for TSC-22. TSC-22 D2 (TSC22 domain family protein 2) also known as TILZ4 is a 780 amino acid protein that is related to TSC-22 and is involved in adaptation of renal cells to hypertonicity, suggesting a possible role in signal transduction. Three isoforms exist due to alternative splicing events.

REFERENCES

- Hino, S., Kawamata, H., Uchida, D., Omotehara, F., Miwa, Y., Begum, N.M., Yoshida, H., Fujimori, T. and Sato, M. 2000. Nuclear translocation of TSC-22 (TGF β -stimulated clone-22) concomitant with apoptosis: TSC-22 as a putative transcriptional regulator. *Biochem. Biophys. Res. Commun.* 278: 659-664.
- Hino, S., Kawamata, H., Omotehara, F., Uchida, D., Begum, N.M., Yoshida, H., Sato, M. and Fujimori, T. 2002. Leucine zipper structure of TSC-22 (TGF β stimulated clone-22) markedly inhibits the anchorage-independent growth of salivary gland cancer cells. *Oncol. Rep.* 9: 371-374.
- Gupta, R.A., Sarraf, P., Brockman, J.A., Shappell, S.B., Raftery, L.A., Willson, T.M. and DuBois, R.N. 2003. Peroxisome proliferator-activated receptor γ and transforming growth factor β pathways inhibit intestinal epithelial cell growth by regulating levels of TSC-22. *J. Biol. Chem.* 278: 7431-7438.
- Uchida, D., Omotehara, F., Nakashiro, K., Tateishi, Y., Hino, S., Begum, N.M., Fujimori, T. and Kawamata, H. 2003. Posttranscriptional regulation of TSC-22 (TGF β -stimulated clone-22) gene by TGF β 1. *Biochem. Biophys. Res. Commun.* 305: 846-854.
- Shostak, K.O., Dmitrenko, V.V., Garifulin, O.M., Rozumenko, V.D., Khomenko, O.V., Zozulya, Y.A., Zehetner, G. and Kavsan, V.M. 2003. Downregulation of putative tumor suppressor gene TSC-22 in human brain tumors. *J. Surg. Oncol.* 82: 57-64.
- Kawamata, H., Fujimori, T. and Imai, Y. 2004. TSC-22 (TGF β stimulated clone-22): a novel molecular target for differentiation-inducing therapy in salivary gland cancer. *Curr. Cancer Drug Targets* 4: 521-529.
- Daouti, S., Latario, B., Nagulapalli, S., Buxton, F., Uziel-Fusi, S., Chirn, G.W., Bodian, D., Song, C., Labow, M., Lotz, M., Quintavalla, J. and Kumar, C. 2005. Development of comprehensive functional genomic screens to identify novel mediators of osteoarthritis. *Osteoarthr. Cartil.* 13: 508-518.
- Shostak, K.O., Dmitrenko, V.V., Vudmaska, M.I., Naidenov, V.G., Beletskii, A.V., Malisheva, T.A., Semenova, V.M., Zozulya, Y.P., Demotes-Mainard, J., Kavsan, V.M. 2005. Patterns of expression of TSC-22 protein in astrocytic gliomas. *Exp. Oncol.* 27: 314-318.

CHROMOSOMAL LOCATION

Genetic locus: TSC22D2 (human) mapping to 3q25.1; Tsc22d2 (mouse) mapping to 3 D.

SOURCE

TSC-22 D2 (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of TSC-22 D2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-55406 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TSC-22 D2 (C-20) is recommended for detection of TSC-22 D2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

TSC-22 D2 (C-20) is also recommended for detection of TSC-22 D2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for TSC-22 D2 siRNA (h): sc-63171, TSC-22 D2 siRNA (m): sc-63172, TSC-22 D2 shRNA Plasmid (h): sc-63171-SH, TSC-22 D2 shRNA Plasmid (m): sc-63172-SH, TSC-22 D2 shRNA (h) Lentiviral Particles: sc-63171-V and TSC-22 D2 shRNA (m) Lentiviral Particles: sc-63172-V.

Molecular Weight of TSC-22 D2: 79 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz MarkerTM compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz MarkerTM Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruzTM Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.