

# NMNAT-3 (L-15): sc-55871

## BACKGROUND

NMNAT proteins are essential cofactors involved in the fundamental processes of cell metabolism. They belong to the eukaryotic NMN adenylyltransferase family. NMNATs participate in the synthesis of NAD<sup>+</sup> by catalyzing the condensation of nicotinamide mononucleotide and ATP. The presence of magnesium and other divalent cations increases their enzymatic activity. The interaction of NMNATs with nuclear proteins is likely to be modulated by phosphorylation. NMNAT proteins contain at least three potential phosphorylation sites and may act as substrates for nuclear kinases. NMNAT-3 (nicotinamide mononucleotide adenylyltransferase-3), also designated PNAT3, is a 252 amino acid protein that localizes to the mitochondria. Highly expressed in the spleen and lungs, NMNAT-3 is able to form homotetramers. Two isoforms exist due to alternative splicing events.

## REFERENCES

1. Sestini, S., et al. 2000. Enzyme activities leading to NAD synthesis in human lymphocytes. *Arch. Biochem. Biophys.* 379: 277-282.
2. Raffaelli, N., et al. 2002. Identification of a novel human nicotinamide mononucleotide adenylyltransferase. *Biochem. Biophys. Res. Commun.* 297: 835-840.
3. Berger, F., et al. 2005. Subcellular compartmentation and differential catalytic properties of the three human nicotinamide mononucleotide adenylyltransferase isoforms. *J. Biol. Chem.* 280: 36334-36341.
4. Online Mendelian Inheritance in Man, OMIM™. 2004. Johns Hopkins University, Baltimore, MD. MIM Number: 608702. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: *Nmnat3* (mouse) mapping to 9 E3.3.

## SOURCE

NMNAT-3 (L-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of NMNAT-3 of rat origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-55871 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

## APPLICATIONS

NMNAT-3 (L-15) is recommended for detection of NMNAT-3 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NMNAT-3 (L-15) is also recommended for detection of NMNAT-3 in additional species, including equine, canine and avian.

Suitable for use as control antibody for NMNAT-3 siRNA (m): sc-62696, NMNAT-3 shRNA Plasmid (m): sc-62696-SH and NMNAT-3 shRNA (m) Lentiviral Particles: sc-62696-V.


Molecular Weight of NMNAT-3: 28 kDa.

Positive Controls: NMNAT-3 (m): 293T Lysate: sc-122083.


## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



NMNAT-3 (L-15): sc-55871. Western blot analysis of NMNAT-3 expression in non-transfected: sc-117752 (A) and mouse transfected: sc-122083 (B) 293T whole cell lysates.



NMNAT-3 (L-15): sc-55871. Western blot analysis of NMNAT-3 expression in non-transfected: sc-117752 (A) and mouse NMNAT-3 transfected: sc-122084 (B) 293T whole cell lysates.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

**MONOS**  
Satisfaction  
Guaranteed

Try **NMNAT-3 (D-10): sc-390433** or **NMNAT-3 (B-9): sc-398848**, our highly recommended monoclonal alternatives to NMNAT-3 (L-15).