

Esophagin (S-19): sc-55941

BACKGROUND

Esophagin, also known as small proline rich protein 3 (SPR3) or cornifin β , belongs to the cornifin family of cornified-envelope structural proteins. It is expressed in mucosal epithelia such as esophagus and tongue and is strongly induced during epidermal keratinocyte differentiation. Due to its highly inducible nature, Esophagin is considered a marker of squamous differentiation. Esophagin serves as a cross-linking protein within the cornified cell envelope and may play a role in the maintenance of normal esophageal epithelial homeostasis. It shares significant homology with the related proteins, SPRR1 and SPRR2. Esophagin is typically not expressed in healthy human epithelium, but its expression is upregulated in numerous hyperproliferative disorders of the skin. In contrast, its expression is dramatically downregulated in esophageal squamous cell carcinoma.

REFERENCES

- Steinert, P.M., et al. 2000. Transglutaminase crosslinking and structural studies of the human small proline rich 3 protein. *Cell Death Differ.* 6: 916-930.
- Smolinski, K.N., et al. 2002. Activation of the Esophagin promoter during esophageal epithelial cell differentiation. *Genomics* 79: 875-880.
- Katou, F., et al. 2003. Differential expression of cornified cell envelope precursors in normal skin, intraorally transplanted skin and normal oral mucosa. *Br. J. Dermatol.* 148: 898-905.
- Kimos, M.C., et al. 2004. Esophagin and proliferating cell nuclear antigen (PCNA) are biomarkers of human esophageal neoplastic progression. *Int. J. Cancer* 111: 415-417.
- Lehr, E., et al. 2004. Infection with human papillomavirus alters expression of the small proline rich proteins 2 and 3. *J. Med. Virol.* 72: 478-483.
- Vitorino, R., et al. 2006. Two-dimensional electrophoresis study of *in vitro* pellicle formation and dental caries susceptibility. *Eur. J. Oral Sci.* 114: 147-153.
- Cai, K., et al. 2006. Decreased polyphenol transport across cultured intestinal cells by a salivary proline-rich protein. *Biochem. Pharmacol.* 71: 1570-1580.
- Stewart, M.G. and Banks, D. 2006. Enhancement of long-term memory retention by colostrin in one-day-old chicks trained on a weak passive avoidance learning paradigm. *Neurobiol. Learn. Mem.* 86: 66-71.
- Zhang, Y., et al. 2007. Exogenous expression of Esophagin/SPRR3 attenuates the tumorigenicity of esophageal squamous cell carcinoma cells via promoting apoptosis. *Int. J. Cancer* 122: 260-266.

CHROMOSOMAL LOCATION

Genetic locus: *Spr3* (mouse) mapping to 3 F1.

SOURCE

Esophagin (S-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Esophagin of mouse origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-55941 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Esophagin (S-19) is recommended for detection of Esophagin of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Esophagin siRNA (m): sc-62283.

Molecular Weight of Esophagin: 37 kDa.

Positive Controls: NIH/3T3 whole cell lysate: sc-2210, 3T3-L1 cell lysate: sc-2243 or F9 cell lysate: sc-2245.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.