



APIP (19F461): sc-56000

BACKGROUND

APIP (Apaf-1-interacting protein), also known as APIP2, is a member of the aldolase class II family and has a highly conserved C-terminal from *C. elegans* to humans. It is ubiquitously expressed, with high expression levels in heart, kidney and skeletal muscle. Alternative splicing produces two isoforms of APIP. Isoform 1 is the full length, 242 amino acid protein; isoform 2 is missing residues 1-38 and contains a distinct sequence for amino acids 39-53. APIP plays an important role in preventing muscle ischemic damage. It suppresses hypoxia-induced cell death by inducing the activation of Akt and ERK 1/2, which are responsible for inhibition of caspase-9 via phosphorylation, and competing with caspase-9 to bind the caspase recruitment domain (CARD) of Apaf-1. Through these mechanisms, APIP negatively regulates the activation of caspase-9 and Apaf-1-mediated cell death.

REFERENCES

1. Lai, C.H., Chou, C.Y., Chang, L.Y., Liu, C.S. and Lin, W. 2000. Identification of novel human genes evolutionarily conserved in *Caenorhabditis elegans* by comparative proteomics. *Genome Res.* 10: 703-713.
2. Cao, G., Xiao, M., Sun, F., Xiao, X., Pei, W., Li, J., Graham, S.H., Simon, R.P. and Chen, J. 2004. Cloning of a novel Apaf-1-interacting protein: a potent suppressor of apoptosis and ischemic neuronal cell death. *J. Neurosci.* 24: 6189-6201.
3. Cho, D.H., Hong, Y.M., Lee, H.J., Woo, H.N., Pyo, J.O., Mak, T.W. and Jung, Y.K. 2004. Induced inhibition of ischemic/hypoxic injury by APIP, a novel Apaf-1-interacting protein. *J. Biol. Chem.* 279: 39942-39950.
4. Anichini, A., Mortarini, R., Sensi, M. and Zanon, M. 2006. Apaf-1 signaling in human melanoma. *Cancer Lett.* 238: 168-179.
5. Schafer, Z.T. and Kornbluth, S. 2006. The apoptosome: physiological, developmental and pathological modes of regulation. *Dev. Cell* 10: 549-561.
6. Cho, D.H., Lee, H.J., Kim, H.J., Hong, S.H., Pyo, J.O., Cho, C. and Jung, Y.K. 2007. Suppression of hypoxic cell death by APIP-induced sustained activation of Akt and ERK 1/2. *Oncogene* 26: 2809-2814.

CHROMOSOMAL LOCATION

Genetic locus: APIP (human) mapping to 11p13.

SOURCE

APIP (19F461) is a mouse monoclonal antibody raised against full length APIP of human origin.

PRODUCT

Each vial contains 100 µg IgG_{2b} in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

APIP (19F461) is recommended for detection of APAF1-interacting protein of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for APIP siRNA (h): sc-61976, APIP shRNA Plasmid (h): sc-61976-SH and APIP shRNA (h) Lentiviral Particles: sc-61976-V.

Molecular Weight of APIP: 27 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.