

HSV-1/2 gB (10B7): sc-56987

BACKGROUND

Membrane fusion is crucial for the entry, spread and formation of enveloped viruses, such as herpes simplex virus, and is mediated by envelope glycoproteins. Two serotypes of the herpes simplex virus, HSV-1 (also known as type 1 or oral) and HSV-2 (type 2 or genital), have been shown to encode at least ten glycoproteins, four of which are necessary and sufficient to facilitate fusion. These four glycoproteins include glycoprotein B (gB), glycoprotein D (gD), glycoprotein H (gH) and glycoprotein L (gL). The fusion event is dependent upon the expression of a gD receptor on target cell membranes and does not require the presence of cell-surface glycosaminoglycans. gB is essential for viral growth as gB free virions are able to bind cells but not to synthesize virus-specific proteins upon infection. HSV-1 and HSV-2 gB exist as a homodimer which may be linked by disulfide bonds. HSV-1 gB is a 904 amino acid protein with an extracellular domain consisting of amino acids 31-730 and a cytoplasmic domain consisting of amino acids 796-904. HSV-2 gB is a protein also containing 904 amino acids with amino acids 23 to 727 making up the extracellular domain and amino acids 793 to 904 making up the cytoplasmic domain.

REFERENCES

1. Cai, W.H., et al. 1988. Role of glycoprotein B of herpes simplex virus type 1 in viral entry and cell fusion. *J. Virol.* 62: 2596-2604.
2. Bystricka, M., et al. 1991. Type-common and type-specific monoclonal antibodies to herpes simplex virus types 1 and 2. *Acta Virol.* 35: 152-164.
3. Slomka, M.J. 1996. Seroepidemiology and control of genital herpes: the value of type-specific antibodies to herpes simplex virus. *Commun. Dis. Rep. CDR Rev.* 6: R41-R45.
4. Bystricka, M., et al. 1997. Monoclonal antibodies to the distinct antigenic sites on glycoproteins C and B and their protective abilities in herpes simplex virus infection. *Acta Virol.* 41: 5-12.
5. Turner, A., et al. 1998. Glycoproteins gB, gD, and gHgL of herpes simplex virus type 1 are necessary and sufficient to mediate membrane fusion in a COS cell transfection system. *J. Virol.* 72: 873-875.
6. Bystricka, M., et al. 1999. Monoclonal antibodies suitable for type-specific identification of herpes simplex viruses by a rapid culture assay. *Acta Virol.* 43: 399-402.

SOURCE

HSV-1/2 gB (10B7) is a mouse monoclonal antibody raised against glycoprotein B of herpesvirus origin.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

HSV-1/2 gB (10B7) is recommended for detection of glycoprotein B of HSV-1 and HSV-2 by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SELECT PRODUCT CITATIONS

1. Chin, L.W., et al. 2010. Anti-herpes simplex virus effects of berberine from *Coptidis rhizoma*, a major component of a Chinese herbal medicine, Ching-Wei-San. *Arch. Virol.* 155: 1933-1941.
2. MacLeod, D.T., et al. 2013. HSV-1 exploits the innate immune scavenger receptor MARCO to enhance epithelial adsorption and infection. *Nat. Commun.* 4: 1963.
3. Heilingloh, C.S., et al. 2015. L particles transmit viral proteins from herpes simplex virus 1-infected mature dendritic cells to uninfected bystander cells, inducing CD83 downmodulation. *J. Virol.* 89: 11046-11055.
4. Petro, C., et al. 2015. Herpes simplex type 2 virus deleted in glycoprotein D protects against vaginal, skin and neural disease. *Elife* 4: e06054.
5. Fabiani, M., et al. 2017. A novel method to titrate herpes simplex virus-1 (HSV-1) using laser-based scanning of near-infrared fluorophores conjugated antibodies. *Front. Microbiol.* 8: 1085.
6. Watanabe, K., et al. 2018. High-sensitivity virus and mycoplasma screening test reveals high prevalence of parvovirus B19 infection in human synovial tissues and bone marrow. *Stem Cell Res. Ther.* 9: 80.
7. Zinser, E., et al. 2018. A new promising candidate to overcome drug resistant herpes simplex virus infections. *Antiviral Res.* 149: 202-210.
8. Michaelis, M., et al. 2019. Omeprazole increases the efficacy of acyclovir against herpes simplex virus type 1 and 2. *Front. Microbiol.* 10: 2790.
9. Grosche, L., et al. 2020. Herpes simplex virus type-2 paralyzes the function of monocyte-derived dendritic cells. *Viruses* 12: 112.
10. Grabowska, K., et al. 2020. αherpesvirus gB homologs are targeted to extracellular vesicles, but they differentially affect MHC class II molecules. *Viruses* 12: 429.
11. Sabatino, M., et al. 2020. Experimental data based machine learning classification models with predictive ability to select *in vitro* active antiviral and non-toxic essential oils. *Molecules* 25: 2452.
12. Birzer, A., et al. 2020. HSV-1 modulates IL-6 receptor expression on human dendritic cells. *Front. Immunol.* 11: 1970.
13. Birzer, A., et al. 2020. Mass spectrometric characterization of HSV-1 L-particles from human dendritic cells and BHK21 cells and analysis of their functional role. *Front. Microbiol.* 11: 1997.
14. Brun, P., et al. 2021. Persistent herpes simplex virus type 1 infection of enteric neurons triggers CD8⁺ T cell response and gastrointestinal neuromuscular dysfunction. *Front. Cell. Infect. Microbiol.* 11: 615350.

RESEARCH USE

For research use only, not for use in diagnostic procedures.