SANTA CRUZ BIOTECHNOLOGY, INC.

EBV Bam HI Z (BDI506): sc-57720



BACKGROUND

Epstein-Barr virus (EBV), also designated human herpesvirus 4 (HHV-4), is a member of the herpesvirus family and is one of the most common human viruses, infecting about 90% of the population. EBV infects only B lymphocytes and, though often asymptomatic, it can cause infectious mononucleosis, a disease characterized by fatigue, fever, sore throat and muscle soreness. The linear genome of EBV circularizes once it enters the cell and exists there as an episome. EBV may play a role in the development of both Burkitt lymphoma, a disease in which a tumor can form on the mandible or maxilla, and nasopharyngeal carcinoma, a tumor found in the upper respiratory tract, most commonly in the nasopharynx. Bam HI E, K, N and Z represent regions in the in EBV genome. Bam HI Z is an immediate early open reading frame that contains the BZLF1 gene, a key gene for EBV replicative cycle. BZLF1 represents a switch from latent infection to lytic EBV infection.

REFERENCES

- Luka, J., Miller, G., Jörnvall, H. and Pearson, G.R. 1986. Characterization of the restricted component of Epstein-Barr virus early antigens as a cytoplasmic filamentous protein. J. Virol. 58: 748-756.
- Goldschmidts, W.L., Ginsburg, M. and Pearson, G.R. 1989. Neutralization of Epstein-Barr virus-induced ribonucleotide reductase with antibody to the major restricted early antigen polypeptide. Virology 170: 330-333.
- Gorgievski-Hrisoho, M., Hinderer, W., Nebel-Schickel, H., Horn, J., Vornhagen, R., Sonneborn, H.H., Wolf, H. and Siegl, G. 1990. Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. J. Clin. Microbiol. 28: 2305-2311.
- Oikawa, O. 1996. Studies on tissue distribution and expression of Epstein-Barr virus using polymerase chain reaction. Hokkaido Igaku Zasshi 70: 729-742.
- Ruf, I.K., Rhyne, P.W., Yang, H., Borza, C.M., Hutt-Fletcher, L.M., Cleveland, J.L. and Sample, J.T. 1999. Epstein-Barr virus regulates c-Myc, apoptosis, and tumorigenicity in Burkitt lymphoma. Mol. Cell. Biol. 19: 1651-1660.
- Gan, Y.Y., Fones-Tan, A., Chan, S.H. and Gan, L.H. 2001. Epstein-Barr viral antigens used in the diagnosis of nasopharyngeal carcinoma. J. Biomed. Sci. 3: 159-169.
- Noorali, S., Pervez, S., Moatter, T., Soomro, I.N., Kazmi, S.U., Nasir, M.I. and Smith, J.L. 2003. Characterization of T cell non-Hodgkin's lymphoma and its association with Epstein-Barr virus in Pakistani patients. Leuk. Lymphoma 44: 807-813.
- 8. Spender, L.C., Lucchesi, W., Bodelon, G., Bilancio, A., Karstegl, C.E., Asano, T., Dittrich-Breiholz, O., Kracht, M., Vanhaesebroeck, B. and Farrell, P.J. 2006. Cell target genes of Epstein-Barr virus transcription factor EBNA-2: induction of the $p55\alpha$ regulatory subunit of Pl3-kinase and its role in survival of EREB2.5 cells. J. Gen. Virol. 87: 2859-2867.

SOURCE

EBV Bam HI Z (BDI506) is a mouse monoclonal antibody raised against EBV Bam HI Z from an infected cell lysate.

PRODUCT

Each vial contains 100 $\mu g~lg G_2$ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

EBV Bam HI Z (BDI506) is recommended for detection of EBV Bam HI Z of Epstein-Barr Virus origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation $[1-2 \mu g \text{ per } 100-500 \mu g \text{ of total protein (1 ml of cell lysate)] and immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).$

Molecular Weight of EBV Bam HI Z: 35 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker[™] compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-mouse IgG-FITC: sc-2010 (dilution range: 1:100-1:400) or goat anti-mouse IgG-TR: sc-2781 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.