

# MRP5 (P-20): sc-5781

## BACKGROUND

The two members of the large family of ABC transporters known to confer multidrug resistance in human cancer cells are the MDR1 P-glycoprotein and the multidrug-resistance protein MRP1. MRP1 is an integral membrane protein that contains an MDR-like core, an N-terminal membrane-bound region and a cytoplasmic linker, and it is expressed in various cerebral cells, as well as in lung, testis and peripheral blood. The MRP gene family also includes MRP2, which is alternatively designated cMOAT (for canalicular multispecific organic anion transporter) and MRP3, which are both conjugate export pumps expressed predominantly in hepatocytes. MRP2 localizes exclusively to the apical membrane and is constitutively expressed at a high level in normal liver cells. Conversely, MRP3 localizes to the basolateral membrane where it also mediates the transport of the organic anion S-(2,4-dinitrophenyl-) glutathione toward the basolateral side of the membrane. MRP3 is normally expressed at comparatively lower levels than MRP2 and increases only when secretion across the apical membrane by MRP2 is impaired. MRP6 protein is highly expressed in liver and kidney, whereas MRP4 and MRP5 are detected in various tissues yet at much lower levels of expression.

## CHROMOSOMAL LOCATION

Genetic locus: ABCC5 (human) mapping to 3q27.1; Abcc5 (mouse) mapping to 16 A3.

## SOURCE

MRP5 (P-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of MRP5 of mouse origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-5781 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## APPLICATIONS

MRP5 (P-20) is recommended for detection of MRP5 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MRP5 siRNA (h): sc-35965, MRP5 siRNA (m): sc-35966, MRP5 shRNA Plasmid (h): sc-35965-SH, MRP5 shRNA Plasmid (m): sc-35966-SH, MRP5 shRNA (h) Lentiviral Particles: sc-35965-V and MRP5 shRNA (m) Lentiviral Particles: sc-35966-V.


Molecular Weight of MRP5: 185 kDa.

Positive Controls: A549 cell lysate: sc-2413, SK-N-SH cell lysate: sc-2410 or HeLa whole cell lysate: sc-2200.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



MRP5 (P-20): sc-5781. Western blot analysis of MRP5 expression in SK-N-SH (A) and HeLa (B) whole cell lysates.

## SELECT PRODUCT CITATIONS

- Thiebaud, N., et al. 2011. Expression and differential localization of xenobiotic transporters in the rat olfactory neuro-epithelium. *Neurosci. Lett.* 505: 180-185.
- Choi, H.J., et al. 2012. Patterns of gene and metabolite define the effects of extracellular osmolality on kidney collecting duct. *J. Proteome Res.* 11: 3816-3828.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **MRP5 (E-10): sc-376965** or **MRP5 (E-4): sc-390797**, our highly recommended monoclonal alternatives to MRP5 (P-20).