



LCMV (M104): sc-57894

BACKGROUND

Lymphocytic choriomeningitis virus (LCMV) is a member of the family *Arenaviridae* that causes a rodent-borne viral infectious nonbacterial/aseptic meningitis, encephalitis or meningoencephalitis. LCMV is a reverse strand single stranded RNA virus that is naturally spread by the common house mouse. Since the virus is somewhat resistant to drying effects, humans can become infected by inhaling infectious aerosolized particles of rodent urine, feces or saliva; by ingesting food contaminated with the virus; by contamination of mucus membranes with infected body fluids; or by directly exposing cuts or other open wounds to virus-infected blood. The virus is usually only dangerous to those who are immunocompromised.

REFERENCES

1. Doherty, P.C., Ceredig, R. and Allan, J.E. 1988. Immunogenetic analysis of cellular interactions governing the recruitment of T lymphocytes and monocytes in lymphocytic choriomeningitis virus-induced immunopathology. *Clin. Immunol. Immunopathol.* 47: 19-26.
2. Salvato, M., Shimomaye, E. and Oldstone, M.B. 1989. The primary structure of the lymphocytic choriomeningitis virus L gene encodes a putative RNA polymerase. *Virology* 169: 377-384.
3. Childs, J.E., Glass, G.E., Korch, G.W., Ksiazek, T.G. and Leduc, J.W. 1992. Lymphocytic choriomeningitis virus infection and house mouse (*Mus musculus*) distribution in urban Baltimore. *Am. J. Trop. Med. Hyg.* 47: 27-34.
4. Roebroek, R.M., Postma, B.H. and Dijkstra, U.J. 1994. Aseptic meningitis caused by the lymphocytic choriomeningitis virus. *Clin. Neurol. Neurosurg.* 96: 178-180.
5. Barton, L.L. and Hyndman, N.J. 2000. Lymphocytic choriomeningitis virus: reemerging central nervous system pathogen. *Pediatrics* 105: E35.
6. Lee, K.J., Novella, I.S., Teng, M.N., Oldstone, M.B. and de La Torre, J.C. 2000. NP and L proteins of lymphocytic choriomeningitis virus (LCMV) are sufficient for efficient transcription and replication of LCMV genomic RNA analogs. *J. Virol.* 74: 3470-3477.
7. Cornu, T.I. and de la Torre, J.C. 2001. RING finger Z protein of lymphocytic choriomeningitis virus (LCMV) inhibits transcription and RNA replication of an LCMV S-segment minigenome. *J. Virol.* 75: 9415-9426.
8. Sánchez, A.B. and de la Torre, J.C. 2005. Genetic and biochemical evidence for an L polymerase of the prototypic arenavirus lymphocytic choriomeningitis virus. *J. Virol.* 79: 7262-7268.
9. Hoey, J. 2005. Lymphocytic choriomeningitis virus. *CMAJ* 173: 1033.

SOURCE

LCMV (M104) is a mouse monoclonal antibody raised against full length LCMV.

STORAGE

For immediate and continuous use, store at 4° C for up to one month. For sporadic use, freeze in working aliquots in order to avoid repeated freeze/thaw cycles. If turbidity is evident upon prolonged storage, clarify solution by centrifugation.

PRODUCT

Each vial contains 1.0 ml culture supernatant containing IgG₁ with < 0.1% sodium azide.

APPLICATIONS

LCMV (M104) is recommended for detection of LCMV infected cells (e.g. Vero, HeLa) of LCMV origin by Western Blotting (starting dilution to be determined by researcher, dilution range 1:10-1:200), immunofluorescence (starting dilution to be determined by researcher, dilution range 1:10-1:200) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

SELECT PRODUCT CITATIONS

1. Sun, Y., Lu, D., Yin, Y., Song, J., Liu, Y., Hao, W., Qi, F., Zhang, G., Zhang, X., Liu, L., Lin, Z., Liang, H., Zhao, X., Jin, Y. and Yin, Y. 2021. PTEN α functions as an immune suppressor and promotes immune resistance in PTEN-mutant cancer. *Nat. Commun.* 12: 5147.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.