

VAMP-1/2 (4H302): sc-58310

BACKGROUND

Syntaxins were originally thought to be docking proteins, but have more recently been categorized as anchoring proteins that anchor themselves to the cytoplasmic surfaces of cellular membranes. Syntaxins have been shown to bind to various proteins involved in exocytosis, including VAMPs (vesicle-associated membrane proteins), NSF (N-ethylmaleimide-sensitive factor), SNAP 25 (synaptosomal-associated protein of 25kDa), SNAPs (soluble NSF attachment proteins) and Synaptotagmin. VAMPs, also designated synaptobrevins, including VAMP-1 and VAMP-2, and Synaptotagmin, a protein that may function as an inhibitor of exocytosis, are vesicular proteins. SNAPs, including α - and γ -SNAP, are cytoplasmic proteins that bind to a membrane receptor complex composed of VAMP, SNAP 25 and Syntaxin. SNAPs mediate the membrane binding of NSF, which is essential for membrane fusion reactions. An additional protein designated synaptophysin may regulate exocytosis by competing with SNAP 25 and syntaxins for VAMP binding.

REFERENCES

- Elferink, L.A., Peterson, M.R. and Scheller, R.H. 1993. A role for Synaptotagmin (p65) in regulated exocytosis. *Cell* 72: 153-159.
- Bennett, M.K., Garcia-Araras, J.E., Elferink, L.A., Peterson, K., Fleming, A.M., Hazuka, C.D. and Scheller, R.H. 1993. The syntaxin family of vesicular transport receptors. *Cell* 74: 863-873.
- Yamaguchi, K. and Akagawa, K. 1994. Exocytosis relating proteins in the nervous system. *Neurosci. Res.* 20: 289-292.
- Hayashi, T., McMahon, H., Yamasaki, S., Binz, T., Hata, Y., Sudhof, T.C. and Niemann, H. 1994. Synaptic vesicle membrane fusion complex: action of clostridial neurotoxins on assembly. *EMBO J.* 13: 5051-5061.
- Edelmann, L., Hanson, P.I., Chapman, E.R. and Jahn, R. 1995. Synaptobrevin binding to synaptophysin: a potential mechanism for controlling the exocytosis fusion machine. *EMBO J.* 14: 224-231.
- McMahon, H.T. and Sudhof, T.C. 1995. Synaptic core complex of synaptobrevin, Syntaxin, and SNAP 25 forms high affinity α -SNAP binding site. *J. Biol. Chem.* 270: 2213-2217.
- Lin, R.C. and Scheller, R.H. 1997. Structural organization of the synaptic exocytosis core complex. *Neuron* 19: 1087-1094.
- Barnard, R.J., Morgan, A. and Burgoyne, R.D. 1997. Stimulation of NSF ATPase activity by α -SNAP is required for SNARE complex disassembly and exocytosis. *J. Cell Biol.* 139: 875-883.

CHROMOSOMAL LOCATION

Genetic locus: VAMP1 (human) mapping to 12p13.31, VAMP2 (human) mapping to 17p13.1; Vamp1 (mouse) mapping to 6 F3, Vamp2 (mouse) mapping to 11 B3.

SOURCE

VAMP-1/2 (4H302) is a mouse monoclonal antibody raised against synaptic vesicle containing fractions immunoprecipitated from brain homogenates of human origin.

PRODUCT

Each vial contains 200 μ g IgM in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.


APPLICATIONS

VAMP-1/2 (4H302) is recommended for detection of VAMP-1 and VAMP-2 of mouse, rat, human, bovine and porcine and *Xenopus laevis* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)].


Molecular Weight of VAMP-1/2: 18 kDa.

Positive Controls: mouse brain extract: sc-2253, VAMP-2 (h2): 293T Lysate: sc-115139 or U-87 MG cell lysate: sc-2411.

DATA



VAMP-1/2 (4H302): sc-58310. Western blot analysis of VAMP-2 expression in non-transfected: sc-117752 (A) and human VAMP-2 transfected: sc-115139 (B) 293T whole cell lysates.



VAMP-1/2 (4H302): sc-58310. Western blot analysis of VAMP-1/2 expression in U-87 MG whole cell lysate.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.