



C3d (030-08): sc-58929

BACKGROUND

The complement component proteins, C2, C3, C4 and C5, are potent anaphylatoxins that are released during complement activation. Binding of these proteins to their respective G protein-coupled receptors, C3aR, C1R and C5aR, induces proinflammatory events, such as cellular degranulation, smooth muscle contraction, arachidonic acid metabolism, cytokine release, leukocyte activation and cellular chemotaxis. C3d is a terminal degradation product of C3 that plays an important role in modulation of the adaptive immune response through the interaction with complement receptor type 2 (CR2). CR2 is important in the switched-isotype, high-affinity and memory humoral immune responses to T-dependent foreign antigens, as well as in the development of the natural antibody repertoire. This pH- and ionic strength-dependent association of C3d with CR2 represents a link between innate and adaptive immunity.

REFERENCES

1. Nagar, B., Jones, R.G., Diefenbach, R.J., Isenman, D.E. and Rini, J.M. 1998. X-ray crystal structure of C3d: a C3 fragment and ligand for CR2. *Science* 280: 1277-1281.
2. Clemenza, L. and Isenman, D.E. 2000. Structure-guided identification of C3d residues essential for its binding to CR2 (CD21). *J. Immunol.* 165: 3839-3848.
3. Szakonyi, G., Guthridge, J.M., Li, D., Young, K., Holers, V.M. and Chen, X.S. 2001. Structure of CR2 in complex with its C3d ligand. *Science* 292: 1725-1728.
4. Gerl, V.B., Bohl, J., Pitz, S., Stoffelns, B., Pfeiffer, N. and Bhakdi, S. 2002. Extensive deposits of complement C3d and C5b-9 in the choriocapillaris of eyes of patients with diabetic retinopathy. *Invest. Ophthalmol. Vis. Sci.* 43: 1104-1108.
5. Morikis, D. and Lambris, J.D. 2004. The electrostatic nature of C3d-CR2 association. *J. Immunol.* 172: 7537-7547.
6. Nakao, M., Miura, C., Itoh, S., Nakahara, M., Okumura, K., Mutsuro, J. and Yano, T. 2004. A complement C3 fragment equivalent to mammalian C3d from the common carp (*Cyprinus carpio*): generation in serum after activation of the alternative pathway and detection of its receptor on the lymphocyte surface. *Fish Shellfish Immunol.* 16: 139-149.
7. Boackle, R.J., Nguyen, Q.L., Leite, R.S., Yang, X. and Vesely, J. 2005. Complement-coated antibody-transfer (CCAT); serum IgA1 antibodies intercept and transport C4 and C3 fragments and preserve IgG₁ deployment (PGD). *Mol. Immunol.* 43: 236-245.
8. Hannan, J.P., Young, K.A., Guthridge, J.M., Asokan, R., Szakonyi, G., Chen, X.S. and Holers, V.M. 2005. Mutational analysis of the CR2 (CR2/CD21)-C3d interaction reveals a putative charged SCR1 binding site for C3d. *J. Mol. Biol.* 346: 845-858.
9. Holers, V.M. 2005. Complement receptors and the shaping of the natural antibody repertoire. *Springer Semin. Immunopathol.* 26: 405-423.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

CHROMOSOMAL LOCATION

Genetic locus: C3 (human) mapping to 19p13.3-p13.2; C3 (mouse) mapping to 17 E1-E3.

SOURCE

C3d (030-08) is a mouse monoclonal antibody raised against full length native C3 of human origin.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

C3d (030-08) is recommended for detection of C3d of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of C3d: 187 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.