

# SDF-1 (C-19): sc-6193

## BACKGROUND

The C-X-C or  $\alpha$  chemokine family is characterized by a pair of cysteine residues separated by a single amino acid and primarily functions as chemoattractants for neutrophils. The C-X-C family includes IL-8, NAP-2, MSGA and stromal cell-derived factor-1, or SDF-1. SDF-1 was originally described as a pre-B cell stimulatory factor, but has now been shown to function as a potent chemoattractant for T cells and monocytes, but not neutrophils. Receptors for the C-X-C family are G protein-coupled, seven-pass, transmembrane domain proteins which include IL-8RA, IL-8RB and fusin (also designated LESTR or CXCR-4). Fusin is highly homologous to the IL-8 receptors, sharing 37% sequence identity at the amino acid level. The IL-8 receptors bind to IL-8, NAP-2 and MSGA, while fusin binds to its cognate ligand, SDF-1. Fusin has been identified as the major co-receptor for T-tropic HIV-1, and SDF-1 has been shown to inhibit HIV-1 infection. Six human SDF-1 isoforms exist due to alternative splicing of CXCL12, the gene encoding SDF-1. Three isoforms are known for mouse and rat.

## CHROMOSOMAL LOCATION

Genetic locus: CXCL12 (human) mapping to 10q11.21; Cxcl12 (mouse) mapping to 6 F1.

## SOURCE

SDF-1 (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of SDF-1 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-6193 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

SDF-1 (C-19) is recommended for detection of all known isoforms of SDF-1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000). SDF-1 (C-19) is also recommended for detection of all known isoforms of SDF-1 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for SDF-1 siRNA (h): sc-39367, SDF-1 siRNA (m): sc-39368, SDF-1 shRNA Plasmid (h): sc-39367-SH, SDF-1 shRNA Plasmid (m): sc-39368-SH, SDF-1 shRNA (h) Lentiviral Particles: sc-39367-V and SDF-1 shRNA (m) Lentiviral Particles: sc-39368-V.

Molecular Weight of SDF-1: 10 kDa.

Positive Controls: RAW 264.7 whole cell lysate: sc-2211.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## SELECT PRODUCT CITATIONS

- Hargreaves, D.C., et al. 2001. A coordinated change in chemokine responsiveness guides plasma cell movements. *J. Exp. Med.* 194: 45-56.
- Callewaere, C., et al. 2007. Cellular and subcellular evidence for neuronal interaction between the chemokine SDF-1/CXCL12 and vasopressin: regulation in the hypothalamo-neurohypophysial system of the Brattleboro rats. *Endocrinology* 149: 310-319.
- Sowinska, A. and Jagodzinski, P.P. 2007. RNA interference-mediated knock-down of DNMT1 and DNMT3B induces CXCL12 expression in MCF-7 breast cancer and AsPC1 pancreatic carcinoma cell lines. *Cancer Lett.* 255: 153-159.
- Borders, A.S., et al. 2007. Macrophage-mediated neuroprotection and neurogenesis in the olfactory epithelium. *Physiol. Genomics* 31: 531-543.
- Hasegawa, T., et al. 2008. Vascular precursors in developing human retina. *Invest. Ophthalmol. Vis. Sci.* 49: 2178-2192.
- Bhattacharyya, B.J., et al. 2008. The chemokine stromal cell-derived factor-1 regulates GABAergic inputs to neural progenitors in the postnatal dentate gyrus. *J. Neurosci.* 28: 6720-6730.
- Kang, H., et al. 2009. Transcript analyses of stromal cell derived factors (SDFs): SDF-2, SDF-4 and SDF-5 reveal a different pattern of expression and prognostic association in human breast cancer. *Int. J. Oncol.* 35: 205-211.
- Kitaori, T., et al. 2009. Stromal cell-derived factor 1/CXCR4 signaling is critical for the recruitment of mesenchymal stem cells to the fracture site during skeletal repair in a mouse model. *Arthritis Rheum.* 60: 813-823.
- Zhu, Y., et al. 2009. SDF1/CXCR4 signalling regulates two distinct processes of precerebellar neuronal migration and its depletion leads to abnormal pontine nuclei formation. *Development* 136: 1919-1928.
- Merino, J., et al. 2009. Splanchnic-brain crosstalk mediated by chemokines in portal hypertensive rats. *Int. J. Gastroenterol.* E-published.
- Martin, S.K., et al. 2010. Hypoxia-inducible factor-2 is a novel regulator of aberrant CXCL12 expression in multiple myeloma plasma cells. *Haematologica* 95: 776-784.
- Xing, D., et al. 2011. Hypoxia and hypoxia-inducible factor in the burn wound. *Wound Repair Regen.* 19: 205-213.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.


 MONOS  
 Satisfaction  
 Guaranteed

Try **SDF-1 (P-159X): sc-74271**, our highly recommended monoclonal alternative to SDF-1 (C-19).