

IL-13R α 1 siRNA (h): sc-63337

BACKGROUND

The Th2 cytokine interleukin-13 (IL-13) plays a critical role in allergen-induced airway hyper-responsiveness (AHR). Two different receptors exist for IL-13, designated IL-13R α 1 and 2. IL-13R α 1 exists as a heterodimer of IL-13R α 1 and IL-4R α as a signaling subunit, whereas IL-13R α 2 acts as a decoy receptor for IL-13. Furthermore, TNF α or IL-4 stimulation induces IL-13R α 2 upregulation, while IL-13R α 1 is constitutively expressed. Cell surface localization of IL-13R α 2 abrogates IL-13 signaling, thus IL-13 induced translocation of the receptor from the cytoplasm provides a mechanism for negative-feedback of IL-13 signaling. The human IL-13R gene maps to chromosome Xq24. IL-13R α 1 expression is predominant in B cells, monocytes and T cells, whereas IL-13R α 2 expression is highest in glioma cells.

REFERENCES

- Guo, J., Apiou, F., Mellerin, M.P., Leneau, B., Jacques, Y. and Minvielle, S. 1997. Chromosome mapping and expression of the human interleukin-13 receptor. *Genomics* 42: 141-145.
- Park, J.W., Taube, C., Yang, E.S., Joetham, A., Balhorn, A., Takeda, K., Miyahara, N., Dakhama, A., Donaldson, D.D. and Gelfand, E.W. 2003. Respiratory syncytial virus-induced airway hyperresponsiveness is independent of IL-13 compared with that induced by allergen. *J. Allergy Clin. Immunol.* 112: 1078-1087.
- Yoshikawa, M., Nakajima, T., Tsukidate, T., Matsumoto, K., Iida, M., Otori, N., Haruna, S., Moriyama, H. and Saito, H. 2003. The negative-feedback regulation of the IL-13 signal by the IL-13 receptor α 2 chain in bronchial epithelial cells. *Cytokine* 24: 293-303.
- Yasunaga, S., Yuyama, N., Arima, K., Tanaka, H., Toda, S., Maeda, M., Matsui, K., Goda, C., Yang, Q., Sugita, Y., Nagai, H. and Izuhara, K. 2003. TNF α and IL-4 regulate expression of IL-13 receptor α 2 on human fibroblasts. *Biochem. Biophys. Res Commun.* 312: 1248-1255.

CHROMOSOMAL LOCATION

Genetic locus: IL13RA1 (human) mapping to Xq24.

PRODUCT

IL-13R α 1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see IL-13R α 1 shRNA Plasmid (h): sc-63337-SH and IL-13R α 1 shRNA (h) Lentiviral Particles: sc-63337-V as alternate gene silencing products.

For independent verification of IL-13R α 1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-63337A, sc-63337B and sc-63337C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

IL-13R α 1 siRNA (h) is recommended for the inhibition of IL-13R α 1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

IL-13R α 1 (D-2): sc-398831 is recommended as a control antibody for monitoring of IL-13R α 1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor IL-13R α 1 gene expression knockdown using RT-PCR Primer: IL-13R α 1 (h)-PR: sc-63337-PR (20 μ l, 484 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.