

# CD20 (B-H20): sc-65298

## BACKGROUND

CD20 is a leukocyte surface antigen consisting of four transmembrane regions and cytoplasmic N- and C-termini. The cytoplasmic domain of CD20 contains multiple phosphorylation sites leading to additional isoforms. CD20 is expressed primarily on B cells but has also been detected on both normal and neoplastic T cells. CD20 functions as a calcium-permeable cation channel, and it is known to accelerate the G<sub>0</sub> to G<sub>1</sub> progression induced by IGF-1. CD20 is activated by the IGF-1 receptor via the  $\alpha$  subunits of the heterotrimeric G proteins. Activation of CD20 significantly increases DNA synthesis and is thought to involve basic helix-loop-helix leucine zipper transcription factors.

## REFERENCES

1. Tedder, T.F., et al. 1994. CD20: a regulator of cell-cycle progression of B lymphocytes. *Immunol. Today* 15: 450-454.
2. Schlossman, S., ed. 1995. *Leucocyte Typing V*. New York: Oxford University Press.
3. Szollosi, J., et al. 1996. Supramolecular complexes of MHC class I, MHC class II, CD20, and tetraspan molecules (CD53, CD81, and CD82) at the surface of a B cell line JY. *J. Immunol.* 157: 2939-2946.
4. Algino, K.M., et al. 1996. CD20 (pan-B cell antigen) expression on bone marrow-derived T cells. *Am. J. Clin. Pathol.* 106: 78-81.
5. Smiers, F.J., et al. 1996. CD20 and CD40 mediated mitogenic responses in B-lineage acute lymphoblastic leukaemia. *Br. J. Haematol.* 93: 125-130.
6. Himmelmann, A., et al. 1997. PU.1/Pip and basic helix loop helix zipper transcription factors interact with binding sites in the CD20 promoter to help confer lineage- and stage-specific expression of CD20 in B lymphocytes. *Blood* 90: 3984-3995.
7. Kanzaki, M., et al. 1997. Activation of the calcium-permeable cation channel CD20 expressed in BALB/c 3T3 cells by Insulin-like growth factor-1. *J. Biol. Chem.* 272: 4964-4969.
8. Mason, D., ed. 2002. *Leucocyte Typing VII*. New York: Oxford University Press.

## CHROMOSOMAL LOCATION

Genetic locus: MS4A1 (human) mapping to 11q12.2.

## SOURCE

CD20 (B-H20) is a mouse monoclonal antibody raised against chronic lymphoid leukemia (CLL) cells of human origin.

## PRODUCT

Each vial contains 100  $\mu$ g IgG<sub>2a</sub> in 1.0 ml of PBS with < 0.1% sodium azide, 0.1% gelatin and 1% BSA.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


## APPLICATIONS

CD20 (B-H20) is recommended for detection of CD20 of human origin by flow cytometry (1  $\mu$ g per 1 x 10<sup>6</sup> cells).

Suitable for use as control antibody for CD20 siRNA (h): sc-29972, CD20 shRNA Plasmid (h): sc-29972-SH and CD20 shRNA (h) Lentiviral Particles: sc-29972-V.

Molecular Weight of CD20 isoforms: 33-37 kDa.

## DATA



CD20 (B-H20): sc-65298. Indirect FCM analysis of human peripheral blood leukocytes stained with CD20 (B-H20), followed by PE-conjugated goat anti-mouse IgG: sc-3738. Quadrant markers were set based on the isotype control, normal mouse IgG<sub>2a</sub>: sc-3878.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.