

HSP 65 (5178): sc-65466

BACKGROUND

The heat shock proteins (HSPs), also referred to as molecular chaperones, comprise a group of highly conserved, abundantly expressed proteins with diverse functions, including the assembly and sequestering of multiprotein complexes, transportation of nascent polypeptide chains across cellular membranes and regulation of protein folding. HSP 60s are a ubiquitous class of HSPs that specifically promote the folding and assembly of cellular polypeptides in an ATP-dependent manner. Mycobacterial HSP 65 is a member of the HSP family that has the ability to facilitate the cross-presentation of an exogenous protein by dendritic cells (DCs) to CD8 T cells. HSP 65 has a high degree of sequence homology with mammalian HSP 60. The presence of HSP 65-specific T cells in blood and in atherosclerotic lesions may be a key factor in initiating the development of atherosclerosis and perpetuating the lesions.

REFERENCES

1. Rambukkana, A., Yong, S. and Das, P.K. 1991. Identification of a novel B cell epitope of restricted specificity on the HSP 65 kDa protein of *Mycobacterium tuberculosis*. *FEMS Microbiol. Immunol.* 3: 39-45.
2. Izaki, S., Goto, Y., Kitamura, K. and Nomaguchi, H. 1994. Antibody to 65 kDa stress protein (HSP 65) of *Mycobacterium leprae* in various inflammatory skin diseases. A preliminary report. *Nippon Rai Gakkai Zasshi* 63: 3-11.
3. Afek, A., George, J., Gilburd, B., Rauova, L., Goldberg, I., Kopolovic, J., Harats, D. and Shoenfeld, Y. 2000. Immunization of low-density lipoprotein receptor deficient (LDL-RD) mice with heat shock protein 65 (HSP 65) promotes early atherosclerosis. *J. Autoimmun.* 14: 115-121.
4. Iltanen, S., Rantala, I., Laippala, P., Holm, K., Partanen, J. and Maki, M. 2000. Expression of HSP 65 in jejunal epithelial cells in patients clinically suspect of coeliac disease. *Autoimmunity* 31: 125-132.
5. Lewthwaite, J.C., Coates, A.R., Tormay, P., Singh, M., Mascagni, P., Poole, S., Roberts, M., Sharp, L. and Henderson, B. 2001. *Mycobacterium tuberculosis* chaperonin 60.1 is a more potent cytokine stimulator than chaperonin 60.2 (HSP 65) and contains a CD14-binding domain. *Infect. Immun.* 69: 7349-7355.
6. Hasan, A., Shinnick, T., Mizushima, Y., van der Zee, R. and Lehner, T. 2002. Defining a T cell epitope within HSP 65 in recurrent aphthous stomatitis. *Clin. Exp. Immunol.* 128: 318-325.
7. Qamra, R. and Mande, S.C. 2004. Crystal structure of the 65 kDa heat shock protein, chaperonin 60.2, of *Mycobacterium tuberculosis*. *J. Bacteriol.* 186: 8105-8113.
8. Kim, H., Kim, S.H., Shim, T.S., Kim, M.N., Bai, G.H., Park, Y.G., Lee, S.H., Chae, G.T., Cha, C.Y., Kook, Y.H. and Kim, B.J. 2005. Differentiation of *Mycobacterium* species by analysis of the heat shock protein 65 gene (HSP 65). *Int. J. Syst. Evol. Microbiol.* 55: 1649-1656.
9. Mia, M.Y., Durai, M., Kim, H.R. and Moudgil, K.D. 2005. Heat shock protein 65-reactive T cells are involved in the pathogenesis of non-antigenic dimethyl dioctadecyl ammonium bromide-induced arthritis. *J. Immunol.* 175: 219-227.

SOURCE

HSP 65 (5178) is a mouse monoclonal antibody raised against purified HSP 65 derivative from *Mycobacterium tuberculosis*.

PRODUCT

Each vial contains 100 µg IgG_{2a} in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

HSP 65 (5178) is recommended for detection of HSP 65 of *Mycobacterium tuberculosis* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000).

Molecular Weight of HSP 65: 65 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.