

HNF-6 (C-20): sc-6559

BACKGROUND

HNF-1 (α and β), HNF-3 (α , β and γ), HNF-4 (α and γ), and HNF-6 compose, in part, a homeoprotein family designated the hepatocyte nuclear factor family. The various HNF-1 isoforms regulate transcription of genes in the liver as well as in other tissues such as kidney, small intestine and thymus. HNF-3 α , HNF-3 β and HNF-3 γ regulate the transcription of numerous hepatocyte genes in adult liver. HNF-3 α and HNF-3 β have also been shown to be involved in gastrulation events such as body axis formation. HNF-4 α and HNF-4 γ have been shown to be important for early embryo development. HNF-4 α is expressed in liver, kidney, pancreas, small intestine, testis and colon; and HNF-4 γ is expressed in each of these tissues except liver. HNF-6 has been shown to bind to the promoter of HNF-3 β , which indicates a potential role of HNF-6 in gut endoderm epithelial cell differentiation. Evidence suggests that HNF-6 may also be a transcriptional activator for at least 22 other hepatocyte-enriched genes, including cytochrome P450 2C13 and α -1 antitrypsin.

REFERENCES

1. Bach, I., et al. 1993. More potent transcriptional activators or a trans-dominant inhibitor of the HNF1 homeoprotein family are generated by alternative RNA processing. *EMBO J.* 12: 4229-4242.
2. Kaestner, K.H., et al. 1994. The HNF-3 gene family of transcription factors in mice: gene structure, cDNA sequence, and mRNA distribution. *Genomics* 20: 377-385.
3. Drewes, T., et al. 1996. Human hepatocyte nuclear factor 4 isoforms are encoded by distinct and differentially expressed genes. *Mol. Cell. Biol.* 16: 925-931.

CHROMOSOMAL LOCATION

Genetic locus: ONECUT1 (human) mapping to 15q21.2; Onecut1 (mouse) mapping to 9 D.

SOURCE

HNF-6 (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of HNF-6 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-6559 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-6559 X, 200 μ g/0.1 ml.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

HNF-6 (C-20) is recommended for detection of HNF-6 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HNF-6 (C-20) is also recommended for detection of HNF-6 in additional species, including canine and porcine.

Suitable for use as control antibody for HNF-6 siRNA (h): sc-37936, HNF-6 siRNA (m): sc-37937, HNF-6 shRNA Plasmid (h): sc-37936-SH, HNF-6 shRNA Plasmid (m): sc-37937-SH, HNF-6 shRNA (h) Lentiviral Particles: sc-37936-V and HNF-6 shRNA (m) Lentiviral Particles: sc-37937-V.

HNF-6 (C-20) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

SELECT PRODUCT CITATIONS

1. Hatzis, P., et al. 2001. Regulatory mechanisms controlling human hepatocyte nuclear factor 4 α gene expression. *Mol. Cell. Biol.* 21: 7320-7330.
2. Rowland Jennifer E., et al. 2005. *In vivo* analysis of growth hormone receptor signaling domains and their associated transcripts. *Mol. Cell. Biol.* 25: 66-77.
3. Kim, I., et al. 2005. Functional characterization of evolutionarily conserved DNA regions in forkhead box F1 gene locus. *J. Biol. Chem.* 280: 37908-37916.
4. Kajiyama, Y., et al. 2006. Characterization of distant enhancers and promoters in the albumin- α -fetoprotein locus during active and silenced expression. *J. Biol. Chem.* 281: 30122-30131.
5. Lehner, F., et al. 2007. The hepatocyte nuclear factor 6 (HNF-6) and FOXA2 are key regulators in colorectal liver metastases. *FASEB J.* 21: 1445-1462.


 MONOS
 Satisfation
 Guaranteed

Try **HNF-6 (G-10): sc-376167** or **HNF-6 (B-1): sc-376308**, our highly recommended monoclonal alternatives to HNF-6 (C-20).