

# 26S Proteasome p48A (12A1): sc-65756

## BACKGROUND

The 26S Proteasome is a large complex involved in the intracellular degradation of proteins in eukaryotes. Ubiquitination by E3 ubiquitin ligases targets proteins for degradation by this complex. The 26S Proteasome plays an important role in the regulation of many biological processes. It is composed of over 30 different subunits and contains at least two copies of each subunit. Assembly of this large complex is ATP-dependent. Due to its size, it is fairly unstable and often disassociates into subcomplexes (including a 20S core and two 19S regulatory complexes). The 26S Proteasome p48A (also known as Rpt3 in yeast and S6 in human) is one of the six ATPase base subunits of the 19S regulatory complex. The 19S regulatory complex recognizes ubiquitinated proteins, removes the ubiquitin chains and translocates the proteins to the 20S core for degradation. The base subunits are involved in protein chaperone activity. The 26S Proteasome p48A is required for the recruitment of the p30 (Rpn12/S14) subunit of the 19S lid.

## REFERENCES

- Kurucz, E., Andó, I., Sümegi, M., Hölzl, H., Kapelari, B., Baumeister, W. and Udvardy, A. 2002. Assembly of the *Drosophila* 26S Proteasome is accompanied by extensive subunit rearrangements. *Biochem. J.* 365: 527-536.
- Szlanka, T., Haracska, L., Kiss, I., Deák, P., Kurucz, E., Andó, I., Virágh, E. and Udvardy, A. 2003. Deletion of proteasomal subunit S5a/Rpn10/p54 causes lethality, multiple mitotic defects and overexpression of proteasomal genes in *Drosophila melanogaster*. *J. Cell Sci.* 116: 1023-1033.
- Ueda, M., Matsui, K., Ishiguro, S., Sano, R., Wada, T., Paponov, I., Palme, K. and Okada, K. 2004. The HALTED ROOT gene encoding the 26S Proteasome subunit RPT2a is essential for the maintenance of *Arabidopsis meristems*. *Development* 131: 2101-2111.
- Babbitt, S.E., Kiss, A., Deffenbaugh, A.E., Chang, Y.H., Bailly, E., Erdjument-Bromage, H., Tempst, P., Buranda, T., Sklar, L.A., Baumler, J., Gogol, E. and Skowyra, D. 2005. ATP hydrolysis-dependent disassembly of the 26S Proteasome is part of the catalytic cycle. *Cell* 121: 553-565.
- Lee, D., Ezhkova, E., Li, B., Pattenden, S.G., Tansey, W.P. and Workman, J.L. 2005. The proteasome regulatory particle alters the SAGA coactivator to enhance its interactions with transcriptional activators. *Cell* 123: 423-436.
- Szutorisz, H., Georgiou, A., Tora, L. and Dillon, N. 2006. The proteasome restricts permissive transcription at tissue-specific gene loci in embryonic stem cells. *Cell* 127: 1375-1388.
- Seong, K.M., Baek, J.H., Yu, M.H. and Kim, J. 2007. Rpn13p and Rpn14p are involved in the recognition of ubiquitinated Gcn4p by the 26S Proteasome. *FEBS Lett.* 581: 2567-2573.
- Vernace, V.A., Arnaud, L., Schmidt-Glenewinkel, T. and Figueiredo-Pereira, M.E. 2007. Aging perturbs 26S Proteasome assembly in *Drosophila melanogaster*. *FASEB J.* 21: 2672-2682.
- Wang, X., Chen, C.F., Baker, P.R., Chen, P.L., Kaiser, P. and Huang, L. 2007. Mass spectrometric characterization of the affinity-purified human 26S Proteasome complex. *Biochemistry* 46: 3553-3565.

## SOURCE

26S Proteasome p48A (12A1) is a mouse monoclonal antibody raised against 26S Proteasome purified from embryos of *Drosophila* origin.

## PRODUCT

Each vial contains 200 µg IgG<sub>2a</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

26S Proteasome p48A (12A1) is recommended for detection of p48A subunit of the 19S regulatory base complex of the 26S Proteasome of *Drosophila melanogaster* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)].


Molecular Weight of 26S Proteasome p48A: 48 kDa.

Positive Controls: *Drosophila* embryo tissue extract.

## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

## DATA



26S Proteasome p48A (12A1): sc-65756. Western blot analysis of 26S Proteasome expression in *Drosophila* embryo (A) tissue extract and purified *Drosophila* 26S Proteasome (B).

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.