

Gemin8 (N-19): sc-66380

BACKGROUND

Gemin8 is a 242 amino acid protein encoded by the mouse gene Gemin8. Gemin8, along with Gemin2-7 and unrip, is a major component of the large multiprotein survival of motor neurons (SMN) complex. The survival of motor neurons (SMN) protein, a product of the disease gene of the common neurodegenerative disease spinal muscular atrophy, is also part of the SMN complex. The SMN complex is a modular composition of proteins with SMN, Gemin8, and Gemin7 in its center. The SMN complex functions as an assembly machine for small nuclear ribonucleoproteins (snRNPs)-the major components of the spliceosome. Gemin8 binds directly to SMN and mediates its interaction with the Gemin6/Gemin7 heterodimer. Importantly, the loss of Gemin6, Gemin7, and unrip interaction with SMN as a result of Gemin8 knock-down affects snRNP assembly by impairing the SMN complex association with Sm proteins but not with snRNAs. The Gemin6/Gemin7 complex binds to Sm proteins and might help organize Sm proteins for formation of Sm rings on snRNA targets.

REFERENCES

- Massenet, S., et al. 2002. The SMN complex is associated with snRNPs throughout their cytoplasmic assembly pathway. *Mol. Cell. Biol.* 22: 6533-6541.
- Shpargel, K.B., et al. 2005. Gemin proteins are required for efficient assembly of Sm-class ribonucleoproteins. *Proc. Natl. Acad. Sci. USA* 102: 17372-17377.
- Carissimi, C., et al. 2006. Gemin8 is a novel component of the survival motor neuron complex and functions in small nuclear ribonucleoprotein assembly. *J. Biol. Chem.* 281: 8126-8134.
- Carissimi, C., et al. 2006. Gemin8 is required for the architecture and function of the survival motor neuron complex. *J. Biol. Chem.* 281: 37009-37016.
- Zhang, H., et al. 2006. Multiprotein complexes of the survival of motor neuron protein SMN with Gemin8 traffic to neuronal processes and growth cones of motor neurons. *J. Neurosci.* 26: 8622-8632.
- Otter, S., et al. 2007. A comprehensive interaction map of the human survival of motor neuron (SMN) complex. *J. Biol. Chem.* 282: 5825-5833.
- Gabanella, F., et al. 2007. Ribonucleoprotein assembly defects correlate with spinal muscular atrophy severity and preferentially affect a subset of spliceosomal snRNPs. *PLoS ONE* 2: e921.
- Donadelli, M., et al. 2007. Synergistic inhibition of pancreatic adenocarcinoma cell growth by trichostatin A and gemcitabine. *Biochim. Biophys. Acta* 1773: 1095-1106.

CHROMOSOMAL LOCATION

Genetic locus: GEMIN8 (human) mapping to Xp22.2; Gemin8 (mouse) mapping to X F5.

SOURCE

Gemin8 (N-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of Gemin8 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-66380 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Gemin8 (N-19) is recommended for detection of Gemin8 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Gemin8 (N-19) is also recommended for detection of Gemin8 in additional species, including canine.

Suitable for use as control antibody for Gemin8 siRNA (h): sc-62372, Gemin8 siRNA (m): sc-62373, Gemin8 shRNA Plasmid (h): sc-62372-SH, Gemin8 shRNA Plasmid (m): sc-62373-SH, Gemin8 shRNA (h) Lentiviral Particles: sc-62372-V and Gemin8 shRNA (m) Lentiviral Particles: sc-62373-V.

Molecular Weight of Gemin8: 32 kDa.

Positive Controls: COLO 320DM cell lysate: sc-2226.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.