

FIRE (L-20): sc-66444

BACKGROUND

FIRE is a 689 amino acid protein that maps to human gene EMR4. FIRE is also known as EGF-like module-containing mucin-like hormone receptor-like 4 and is a member of the G protein-coupled receptor 2 family, LN-TM7 sub-family. The epidermal growth factor-transmembrane seven (EGF-TM7) family are proteins that express EGF-like domains at their extracellular N-terminus coupled to a classical seven-transmembrane (TM7) cassette. FIRE is believed to mediate the cellular interaction between myeloid cells and B-cells. FIRE is a membrane bound, multi-pass protein that is proteolytically cleaved into two subunits, an extracellular α subunit and a seven-transmembrane subunit. FIRE contains one GPS (G protein-coupled receptor proteolytic cleavage site) domain and two EGF (epidermal growth factor)-like domains, the second of which mediates the interaction with the putative ligand.

REFERENCES

1. Stacey, M., Chang, G.W., Sanos, S.L., Chittenden, L.R., Stubbs, L., Gordon, S. and Lin, H.H. 2002. EMR4, a novel epidermal growth factor (EGF)-TM7 molecule up-regulated in activated mouse macrophages, binds to a putative cellular ligand on B lymphoma cell line A20. *J. Biol. Chem.* 277: 29283-29293.
2. Rautava, J., Soukka, T., Heikinheimo, K., Miettinen, P.J., Happonen, R.P. and Jaakkola, P. 2003. Different mechanisms of Syndecan-1 activation through a fibroblast-growth-factor-inducible response element (FIRE) in mucosal and cutaneous wounds. *J. Dent. Res.* 82: 382-387.
3. Hamann, J., Kwakkenbos, M.J., de Jong, E.C., Heus, H., Olsen, A.S. and van Lier, R.A. 2003. Inactivation of the EGF-TM7 receptor EMR4 after the pan-homo divergence. *Eur. J. Immunol.* 33: 1365-1371.
4. Leemans, J.C., te Velde, A.A., Florquin, S., Bennink, R.J., de Bruin, K., van Lier, R.A., van der Poll, T. and Hamann, J. 2004. The epidermal growth factor-seven transmembrane (EGF-TM7) receptor CD97 is required for neutrophil migration and host defense. *J. Immunol.* 172: 1125-1131.
5. Kwakkenbos, M.J., Kop, E.N., Stacey, M., Matmati, M., Gordon, S., Lin, H.H. and Hamann, J. 2004. The EGF-TM7 family: a postgenomic view. *Immunogenetics* 55: 655-666.
6. Hamann, J. 2004. The EGF-TM7 family of the rat. *Immunogenetics* 56: 679-681.
7. Caminschi, I., Vandenabeele, S., Sofi, M., McKnight, A.J., Ward, N., Brodnicki, T.C., Toy, T., Lahoud, M., Maraskovsky, E., Shortman, K. and Wright, M.D. 2006. Gene structure and transcript analysis of the human and mouse EGF-TM7 molecule, FIRE. *DNA Seq.* 17: 8-14.

SOURCE

FIRE (L-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of FIRE of rat origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-66444 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

FIRE (L-20) is recommended for detection of FIRE of rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of FIRE: 77 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.