



## ORP-2 (I-16): sc-66569

### BACKGROUND

Members of the oxysterol-binding protein (OSBP) family function as intracellular lipid receptors. OSBPs are involved in lipid metabolism and signal transduction, as well as vesicle transport, and can translocate to the periphery of Golgi membranes when they are bound to oxysterols. ORPs (OSBP-related proteins) belong to a subfamily of OSBPs and consist of ORP-1 and ORP-2. The ORPs have a highly conserved OSBP-type sterol-binding region and a Pleckstrin homology domain. They strongly bind to phosphatidic acid and weakly bind to phosphatidylinositol 3-phosphate. Two isoforms of ORP-2 are produced due to alternative splicing. ORP-2 is expressed only in retina, retinal pigment epithelium choroid, pineal gland and cultured retinal pigment epithelial cells. Overexpression of ORP-2 decreases cell growth and blocks Golgi-derived vesicle transport.

### REFERENCES

- Laitinen, S., Oikkonen, V.M., Ehnholm, C. and Ikonen, E. 1999. Family of human oxysterol binding protein (OSBP) homologues. A novel member implicated in brain sterol metabolism. *J. Lipid Res.* 40: 2204-2211.
- Xu, Y., Liu, Y., Ridgway, N.D. and McMaster, C.R. 2001. Novel members of the human oxysterol-binding protein family bind phospholipids and regulate vesicle transport. *J. Biol. Chem.* 276: 18407-18414.
- Lehto, M., Laitinen, S., Chinetti, G., Johansson, M., Ehnholm, C., Staels, B., Ikonen, E. and Oikkonen, V.M. 2001. The OSBP-related protein family in humans. *J. Lipid Res.* 42: 1203-1213.
- Jaworski, C.J., Moreira, E., Li, A., Lee, R. and Rodriguez, I.R. 2001. A family of human genes containing oxysterol-binding domains. *Genomics* 78: 185-196.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 606731. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Fairn, G.D. and McMaster, C.R. 2005. Identification and assessment of the role of a nominal phospholipid binding region of ORP(oxysterol-binding-protein-related protein 1 short) in the regulation of vesicular transport. *Biochem. J.* 387: 889-896.
- Perry, R.J. and Ridgway, N.D. 2006. Oxysterol-binding protein and vesicle-associated membrane protein-associated protein are required for sterol-dependent activation of the ceramide transport protein. *Mol. Biol. Cell.* 17: 2604-2616.
- Oikkonen, V.M., Johansson, M., Suchanek, M., Yan, D., Hynynen, R., Ehnholm, C., Jauhainen, M., Thiele, C. and Lehto, M. 2006. The OSBP-related proteins (ORPs): global sterol sensors for co-ordination of cellular lipid metabolism, membrane trafficking and signalling processes? *Biochem. Soc. Trans.* 34: 389-391.
- SWISS-PROT/TrEMBL (Q9H1P3). World Wide Web URL: <http://www.expasy.ch/sprot/sprot-top.html>

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### CHROMOSOMAL LOCATION

Genetic locus: OSBPL2 (human) mapping to 20q13.3.

### SOURCE

ORP-2 (I-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of ORP-2 of human origin.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-66569 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

ORP-2 (I-16) is recommended for detection of Oxysterol-binding protein-related protein 2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for ORP-2 siRNA (h): sc-62717.

Molecular Weight of ORP-2: 55 kDa.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.