

FMNL1 (P-17): sc-66760

BACKGROUND

Formin-like protein 1 (FMNL1, formin-related protein, Frl) is a 1,094 amino acid protein encoded by the mouse gene *Fmnl1*. FMNL1 belongs to the formin homology family and has one DAD (diaphanous autoregulatory domain), one FH2 (formin homology 2) domain and one GBD/FH3 (Rho GTPase-binding/formin homology 3) domain. Formins are a conserved class of proteins expressed in all eukaryotes, with known roles in generating cellular actin-based structures. FMNL1 is believed to play a role in the control of cell motility and survival of macrophages. FMNL1 has been found to interact with RAC1, PFN1 and PFN2 and can block apoptotic cell death and inhibit cell adhesion and migration. FMNL1 is located in the cytoplasm and is highly expressed in the spleen, lymph node and bone marrow cells.

REFERENCES

1. Yayoshi-Yamamoto, S., et al. 2000. FRL, a novel formin-related protein, binds to Rac and regulates cell motility and survival of macrophages. *Mol. Cell. Biol.* 20: 6872-6881.
2. Katoh, M. and Katoh, M. 2003. Identification and characterization of human FMNL1, FMNL2 and FMNL3 genes in silico. *Int. J. Oncol.* 22: 1161-1168.
3. Katoh, M. and Katoh, M. 2004. Identification and characterization of the human FMN1 gene in silico. *Int. J. Mol. Med.* 14: 121-126.
4. Harris, E.S., et al. 2004. The mouse formin, FRL α , slows actin filament barbed end elongation, competes with capping protein, accelerates polymerization from monomers and severs filaments. *J. Biol. Chem.* 279: 20076-20087.
5. Favaro, P.M., et al. 2006. High expression of FMNL1 protein in T non-Hodgkin's lymphomas. *Leuk. Res.* 30: 735-738.
6. Gomez, T.S., et al. 2007. Formins regulate the Actin-related protein 2/3 complex-independent polarization of the centrosome to the immunological synapse. *Immunity* 26: 177-190.

CHROMOSOMAL LOCATION

Genetic locus: FMNL1 (human) mapping to 17q21.31.

SOURCE

FMNL1 (P-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of FMNL1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-66760 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4 $^{\circ}$ C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

FMNL1 (P-17) is recommended for detection of FMNL1, also designated Formin-like protein 1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for FMNL1 siRNA (h): sc-62325, FMNL1 shRNA Plasmid (h): sc-62325-SH and FMNL1 shRNA (h) Lentiviral Particles: sc-62325-V.


Molecular Weight of FMNL1: 160 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204, HuT 78 whole cell lysate: sc-2208 or FMNL1 (h): 293T lysate: sc-113684.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



FMNL1 (P-17): sc-66760. Western blot analysis of FMNL1 expression in Jurkat (A), NK-92 (B) and HUT 78 (C) whole cell lysates.



FMNL1 (P-17): sc-66760. Western blot analysis of FMNL1 expression in non-transfected: sc-117752 (A) and human FMNL1 transfected: sc-113684 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.


 MONOS
Satisfaction
Guaranteed

Try **FMNL1 (C-5): sc-390023** or **FMNL1 (A-4): sc-390466**, our highly recommended monoclonal alternatives to FMNL1 (P-17).