

PDZ-RhoGEF (R-295): sc-67024

BACKGROUND

The multidomain (RGS)-containing RhoGEFs represent a family of guanine nucleotide exchange factors that stabilize the nucleotide-free state of small GTPases through their DH/PH domains, leading to the exchange of GDP to GTP. Uniquely, PDZ-RhoGEF, also known as Rho guanine nucleotide exchange factor 11 and ARHGEF11, binds tightly to both nucleotide-free and activated Rho A, therefore playing a role as a primary regulator of Rho A. Mutations within the carboxylate-binding loop of PDZ-RhoGEF result in changes in cell morphology and Actin organization which is likely due to its interaction with MAP-1A (MAP1 light chain LC2). PDZ-RhoGEF also plays a role in B plexin-mediated activation of Rho/Rho kinase signaling, which is implicated in the regulation of axon guidance and cell migration.

REFERENCES

1. Fukuhara, S., et al. 1999. A novel PDZ domain containing guanine nucleotide exchange factor links heterotrimeric G proteins to Rho. *J. Biol. Chem.* 274: 5868-5879.
2. Rumenapp, U., et al. 1999. Rho-specific binding and guanine nucleotide exchange catalysis by KIAA0380, a DBL family member. *FEBS. Lett.* 459: 313-318.
3. Garrard, S.M., et al. 2001. Expression, purification, and crystallization of the RGS-like domain from the Rho nucleotide exchange factor, PDZ-RhoGEF, using the surface entropy reduction approach. *Protein. Expr. Purif.* 21: 412-416.
4. Driessens, M.H., et al. 2002. B plexins activate Rho through PDZ-RhoGEF. *FEBS. Lett.* 529: 168-172.
5. Oleksy, A., et al. 2004. Preliminary crystallographic analysis of the complex of the human GTPase Rho A with the DH/PH tandem of PDZ-RhoGEF. *Acta. Crystallogr. D Biol. Crystallogr.* 60: 740-742.
6. Tanabe, S., et al. 2004. Regulation of RGS-RhoGEFs by G $_{\alpha 12}$ and G $_{\alpha 13}$ proteins. *Methods Enzymol.* 390: 285-294.
7. Longhurst, D.M., et al. 2006. Interaction of PDZ-RhoGEF with microtubule-associated protein 1 light chains: link between microtubules, Actin cytoskeleton, and neuronal polarity. *J. Biol. Chem.* 281: 12030-12040.
8. Online Mendelian Inheritance in Man, OMIM™. 2009. Johns Hopkins University, Baltimore, MD. MIM Number: 605708. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
9. Chen, Z., et al. 2010. Activated Rho A binds to the PH domain of PDZ-RhoGEF: A potential site for autoregulation. *J. Biol. Chem.* 285: 21070-21081.

CHROMOSOMAL LOCATION

Genetic locus: Arhgef11 (mouse) mapping to 3 F1.

SOURCE

PDZ-RhoGEF (R-295) is a rabbit polyclonal antibody raised against amino acids 1228-1522 mapping at the C-terminus of PDZ-RhoGEF of rat origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

PDZ-RhoGEF (R-295) is recommended for detection of PDZ-RhoGEF of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for PDZ-RhoGEF siRNA (m): sc-45824, PDZ-RhoGEF shRNA Plasmid (m): sc-45824-SH and PDZ-RhoGEF shRNA (m) Lentiviral Particles: sc-45824-V.

Molecular Weight of PDZ-RhoGEF: 183 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfation
Guaranteed

Try **PDZ-RhoGEF (20): sc-136469** or **PDZ-RhoGEF (D-9): sc-166740**, our highly recommended monoclonal alternatives to PDZ-RhoGEF (R-295).