

SUZ12 (H-300): sc-67105

BACKGROUND

The Polycomb group (PcG) genes contribute to the maintenance of cell identity, cell cycle regulation and oncogenesis. The mammalian PcG proteins are regulatory proteins important for Hox gene expression, axial skeleton development and the control of proliferation and survival of hematopoietic cells. By inducing changes in chromatin structure, the PcG proteins are part of a cellular memory system that is responsible for gene activity being inherited to progeny cells. PcG proteins silence gene expression through the formation of multimeric protein complexes with different compositions. Manipulating the expression-levels of various PcG proteins in mammalian cell lines results in cellular transformation, which may be a link between the chromatin-associated PcG proteins and cancer. Polycomb protein SUZ12, also designated ChET 9 protein, JJAZ1 or joined to JAZF1 protein, is a nuclear protein belonging to the VEFs (VRN2-EMF2-FIS2-SUZ12) family. The gene maps to chromosome 17q11.2. SUZ12 has been detected at the breakpoints of a certain recurrent chromosomal translocation which has been reported in endometrial stromal sarcoma. It is a component of the PRC2 complex, composed of EED, EZH2, SUZ12, RBBP4 and RBBP7.

REFERENCES

- Weinmann, A.S., et al. 2001. Use of chromatin immunoprecipitation to clone novel E2F target promoters. *Mol. Cell. Biol.* 21: 6820-6832.
- Koontz, J.I., et al. 2001. Frequent fusion of the JAZF1 and JJAZ1 genes in endometrial stromal tumors. *Proc. Natl. Acad. Sci. USA* 98: 6348-6353.
- Cao, R., et al. 2002. Role of Histone H3 Lysine 27 methylation in Polycomb-group silencing. *Science* 298: 1039-1043.
- Kuzmichev, A., et al. 2002. Histone methyltransferase activity associated with a human multiprotein complex containing the enhancer of zeste protein. *Genes Dev.* 16: 2893-2905.
- Kirmizis, A., et al. 2004. Silencing of human polycomb target genes is associated with methylation of Histone H3 Lys 27. *Genes Dev.* 18: 1592-1605.
- Cao, R., et al. 2004. SUZ12 is required for both the histone methyltransferase activity and the silencing function of the EED-EZH2 complex. *Mol. Cell* 15: 57-67.
- Pasini, D., et al. 2004. SUZ12 is essential for mouse development and for EZH2 histone methyltransferase activity. *EMBO J.* 23: 4061-4071.

CHROMOSOMAL LOCATION

Genetic locus: SUZ12 (human) mapping to 17q11.2; Suz12 (mouse) mapping to 11 B5.

SOURCE

SUZ12 (H-300) is a rabbit polyclonal antibody raised against amino acids 440-739 mapping at the C-terminus of SUZ12 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

SUZ12 (H-300) is recommended for detection of SUZ12 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1–2 µg per 100–500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


SUZ12 (H-300) is also recommended for detection of SUZ12 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for SUZ12 siRNA (h): sc-45597, SUZ12 siRNA (m): sc-45598, SUZ12 shRNA Plasmid (h): sc-45597-SH, SUZ12 shRNA Plasmid (m): sc-45598-SH, SUZ12 shRNA (h) Lentiviral Particles: sc-45597-V and SUZ12 shRNA (m) Lentiviral Particles: sc-45598-V.

Molecular Weight of SUZ12: 95 kDa.

Positive Controls: SUZ12 (h): 293 Lysate: sc-111839, HeLa whole cell lysate: sc-2200 or SW480 cell lysate: sc-2219.

DATA



SUZ12 (H-300): sc-67105. Western blot analysis of SUZ12 expression in non-transfected 293: sc-110760 (A), human SUZ12 transfected 293: sc-111839 (B), HeLa (C), SW480 (D) and JEG-3 (E) whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **SUZ12 (D-10): sc-271325**, our highly recommended monoclonal alternative to SUZ12 (H-300).