

μ-protocadherin (M-303): sc-67367

BACKGROUND

The mucin-like protocadherin, μ-protocadherin, is a developmentally regulated, single pass type I transmembrane protein that belongs to the cadherin superfamily. It contains four cadherin-like ectodomains, a triply repeating mucin domain, four SH3 binding regions, N- and O-glycosylation sites and a possible C-terminal PDZ binding domain. μ-protocadherin is expressed in various epithelial tissues and localizes to the apical surface along the brush border of the proximal convoluted tubule. It acts as a calcium-dependent cell adhesion molecule mediating cell aggregation and may play a role in the activation of signaling events. Due to alternative splicing at least four isoforms exist for μ-protocadherin. These isoforms vary in the region containing the mucin-like domains. Only the longest isoform contains the triply repeating mucin domain.

REFERENCES

- Goldberg, M., et al. 2000. μ-protocadherin, a novel developmentally regulated protocadherin with mucin-like domains. *J. Biol. Chem.* 275: 24622-24629.
- Paris, M.J. and Williams, B.R. 2001. Characterization of a 500-kb contig spanning the region between c-Ha-Ras and MUC2 on chromosome 11p15.5. *Genomics* 69: 196-202.
- Goldberg, M., et al. 2002. Identification and expression analysis of the human μ-protocadherin gene in fetal and adult kidneys. *Am. J. Physiol. Renal Physiol.* 283: F454-F463.
- Goldberg, M., et al. 2003. Biallelic expression of HRAS and MUCDHL in human and mouse. *Hum. Genet.* 112: 334-342.
- Wang, Y., et al. 2004. Gene expression profiles and molecular markers to predict recurrence of Dukes' B colon cancer. *J. Clin. Oncol.* 22: 1564-1571.
- Moulton, D.E., et al. 2004. Expression of a novel cadherin in the mouse and human intestine. *Pediatr. Res.* 55: 927-934.
- Redies, C., et al. 2005. δ-protocadherins: unique structures and functions. *Cell. Mol. Life Sci.* 62: 2840-2852.
- Kostadinov, R., et al. 2006. GRSDb: a database of quadruplex forming G-rich sequences in alternatively processed mammalian pre-mRNA sequences. *Nucleic Acids Res.* 34: D119-D124.
- Garrison, W.D., et al. 2006. Hepatocyte nuclear factor 4alpha is essential for embryonic development of the mouse colon. *Gastroenterology* 130: 1207-1220.

CHROMOSOMAL LOCATION

Genetic locus: Mucdhl (mouse) mapping to 7 F5.

SOURCE

μ-protocadherin (M-303) is a rabbit polyclonal antibody raised against amino acids 209-511 mapping within an extracellular domain of μ-protocadherin of mouse origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 μg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

μ-protocadherin (M-303) is recommended for detection of μ-protocadherin of mouse and, to a lesser extent, rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1–2 μg per 100–500 μg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for μ-protocadherin siRNA (h): sc-72286, μ-protocadherin shRNA Plasmid (h): sc-72286-SH and μ-protocadherin shRNA (h) Lentiviral Particles: sc-72286-V.


Molecular Weight of μ-protocadherin isoforms: 110-220 kDa.

Positive Controls: μ-protocadherin (m): 293T Lysate: sc-127848.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



μ-protocadherin (M-3): sc-67367. Western blot analysis of μ-protocadherin expression in non-transfected: sc-117752 (A) and mouse μ-protocadherin transfected: sc-127848 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.