

MFRP (H-194): sc-67381

BACKGROUND

MFRP (membrane frizzled-related protein) is a single pass type II membrane protein with two cubilin (CUB) domains, two LDL-receptor class A domains, a cysteine-rich domain (CRD) and seven N-glycosylation sites. The C-terminal CRD is related to the Wnt-binding domain of the frizzled family of transmembrane proteins, suggesting that MFRP may act as a Wnt receptor. MFRP is predominantly expressed in retinal pigment epithelial cells, co-localizing to the plasma membrane with CTRP5. MFRP interacts with CTRP5 via its CUB domain and may play a role in the development of the eye. Functional MFRP is necessary for photoreceptor maintenance and appears to regulate the axial length of the eye. Mutations in the gene encoding MFRP can affect various structures in the eye and may result in nanophthalmos 2 (NNO2), an eye disorder characterized by extreme hyperopia.

REFERENCES

1. Katoh, M. 2001. Molecular cloning and characterization of MFRP, a novel gene encoding a membrane-type frizzled-related protein. *Biochem. Biophys. Res. Commun.* 282: 116-123.
2. Kameya, S., et al. 2002. Mfrp, a gene encoding a frizzled-related protein, is mutated in the mouse retinal degeneration 6. *Hum. Mol. Genet.* 11: 1879-1886.
3. Jones, S.E. and Jomary, C. 2002. Secreted frizzled-related proteins: searching for relationships and patterns. *Bioessays* 24: 811-820.
4. Pauer, G.J., et al. 2005. Mutation screen of the membrane-type frizzled-related protein (MFRP) gene in patients with inherited retinal degenerations. *Ophthalmic Genet.* 26: 157-161.
5. Sundin, O.H., et al. 2005. Extreme hyperopia is the result of null mutations in MFRP, which encodes a frizzled-related protein. *Proc. Natl. Acad. Sci. USA* 102: 9553-9558.
6. Katoh, M. and Katoh, M. 2005. Identification and characterization of rat Ror1 and Ror2 genes in silico. *Int. J. Mol. Med.* 15: 533-538.
7. Ayala-Ramirez, R., et al. 2006. A new autosomal recessive syndrome consisting of posterior microphthalmos, retinitis pigmentosa, foveoschisis, and optic disc drusen is caused by a MFRP gene mutation. *Mol. Vis.* 12: 1483-1489.
8. Mandal, M.N., et al. 2006. Spatial and temporal expression of MFRP and its interaction with CTRP5. *Invest. Ophthalmol. Vis. Sci.* 47: 5514-5521.
9. Shu, X., et al. 2006. Disease mechanisms in late-onset retinal macular degeneration associated with mutation in C10TNF5. *Hum. Mol. Genet.* 15: 1680-1689.

CHROMOSOMAL LOCATION

Genetic locus: MFRP (human) mapping to 11q23.3; Mfrp (mouse) mapping to 9 A5.1.

SOURCE

MFRP (H-194) is a rabbit polyclonal antibody raised against amino acids 386-579 mapping within a C-terminal extracellular domain of MFRP of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

MFRP (H-194) is recommended for detection of MFRP of human and, to a lesser extent, mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


MFRP (H-194) is also recommended for detection of MFRP in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for MFRP siRNA (h): sc-72373, MFRP siRNA (m): sc-77398, MFRP shRNA Plasmid (h): sc-72373-SH, MFRP shRNA Plasmid (m): sc-77398-SH, MFRP shRNA (h) Lentiviral Particles: sc-72373-V and MFRP shRNA (m) Lentiviral Particles: sc-77398-V.


Molecular Weight of MFRP: 62 kDa.

Positive Controls: TT whole cell lysate: sc-364195.

DATA



MFRP (H-194): sc-67381. Western blot analysis of MFRP expression in non-transfected: sc-110760 (A) and mouse MFRP transfected: sc-178932 (B) 293 whole cell lysates.



MFRP (H-194): sc-67381. Western blot analysis of MFRP expression in TT whole cell lysate.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.