

# Otopiralin (C-18): sc-67728

## BACKGROUND

Otopiralin is an 89 amino acid inner ear-specific protein encoded by the OTOS gene. Otopiralin is synthesized by fibrocytes of spiral limbus and spiral ligament in the cochlea. Fibrocytes are responsible for maintaining inner ear homeostasis and impairment or alteration of these cells may lead to deterioration of auditory function. Degeneration of fibrocytes due to the absence of Otopiralin leads to irreversible deafness in guinea pigs and moderate deafness in mice. Loss of function in hair cells of the inner ear may also be caused by the downregulation of Otopiralin. Otopiralin is conserved from fish to mammals. It shares homology with Gag p30 core shell and SARS of type C retroviruses. One isoform is produced due to alternative splicing.

## REFERENCES

1. Gratton, M.A., Schulte, B.A. and Hazen-Martin, D.J. 1996. Characterization and development of an inner ear type I fibrocyte cell culture. *Hear. Res.* 99: 71-78.
2. Delprat, B., Boulanger, A., Wang, J., Beaudoin, V., Guitton, M.J., Ventéo, S., Dechesne, C.J., Pujol, R., Lavigne-Rebillard, M., Puel, J.L. and Hamel, C.P. 2002. Downregulation of Otopiralin, a novel inner ear protein, causes hair cell degeneration and deafness. *J. Neurosci.* 22: 1718-1725.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607877. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Lavigne-Rebillard, M., Delprat, B., Surget, M.O., Griffoin, J.M., Weil, D., Arbones, M., Vincent, R. and Hamel, C.P. 2003. Gene structure, chromosomal localization and mutation screening of the human gene for the inner ear protein Otopiralin. *Neurogenetics* 4: 137-140.
5. Pompeia, C., Hurle, B., Belyantseva, I.A., Noben-Trauth, K., Beisel, K., Gao, J., Buchoff, P., Wistow, G. and Kachar, B. 2004. Gene expression profile of the mouse organ of Corti at the onset of hearing. *Genomics* 83: 1000-1011.
6. Caravelli, A., Pianese, L., Saulino, C., Di Leva, F., Sequino, L., Cocozza, S., Marciano, E. and Franzé, A. 2004. Downregulation of Otopiralin mRNA in response to acoustic stress in guinea pig. *Hear. Res.* 198: 36-40.
7. Delprat, B., Ruel, J., Guitton, M.J., Hamard, G., Lenoir, M., Pujol, R., Puel, J.L., Brabet, P. and Hamel, C.P. 2005. Deafness and cochlear fibrocyte alterations in mice deficient for the inner ear protein Otopiralin. *Mol. Cell. Biol.* 25: 847-853.

## CHROMOSOMAL LOCATION

Genetic locus: OTOS (human) mapping to 2q37.3; Otos (mouse) mapping to 1 D.

## SOURCE

Otopiralin (C-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Otopiralin of human origin.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-67728 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Otopiralin (C-18) is recommended for detection of Otopiralin of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Otopiralin (C-18) is also recommended for detection of Otopiralin in additional species, including canine.

Suitable for use as control antibody for Otopiralin siRNA (h): sc-62725, Otopiralin siRNA (m): sc-62726, Otopiralin shRNA Plasmid (h): sc-62725-SH, Otopiralin shRNA Plasmid (m): sc-62726-SH, Otopiralin shRNA (h) Lentiviral Particles: sc-62725-V and Otopiralin shRNA (m) Lentiviral Particles: sc-62726-V.

Molecular Weight of Otopiralin: 6 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **Otopiralin (A-12): sc-393239**, our highly recommended monoclonal alternative to Otopiralin (C-18).