SANTA CRUZ BIOTECHNOLOGY, INC.

Cdc24 (yN-18): sc-6794



BACKGROUND

One of the best characterized MAP kinase pathways is the yeast pheromone response pathway. Extracellular pheromones bind to cell surface receptors and lead to the activation of the MAPKKK Ste11 and the subsequent members of this cascade, Ste7, FusS3 (also designated Dac2) and Kss1. These MAP kinases activate Ste12 and Far1, which effect transcriptional and morphological changes necessary for mating. Cdc42, a small GTP-binding protein, is thought to play a critical role in this pathway. Activity of Cdc42 is regulated in part by the guanylyl-nucleotide-exchange factor Cdc24 (also designated Cls4). Cells lacking Cdc24 are unable to respond to pheromone and also show defects in bud-site selection and bud emergence.

REFERENCES

- Errede, B. and Ammerer, G. 1989. STE12, a protein involved in cell-typespecific transcription and signal transduction in yeast, is part of protein-DNA complexes. Genes and Dev. 3: 1349-1361.
- Peter, M., Gartner, A., Horecka, J., Ammerer, G., and Herskowitz, I. 1993. FAR1 links the signal transduction pathway to the cell cycle machinery in yeast. Cell 73: 747-760.
- Zheng, Y., Cerione, R., and Bender, A. 1994. Control of the yeast bud-site assembly GTPase Cdc42. J. Biol. Chem. 269: 2369-2372.
- Ferguson, B., Horecka, J., Printen, J., Schultz, J., Stevenson, B.J., and Sprague, G.F. Jr. 1994. The yeast pheromone response pathway: new insights into signal transmission. Cell. Mol. Biol. Res. 40: 223-228.
- 5. Herskowitz, I. 1995. MAP kinase pathways in yeast: for mating and more. Cell 80: 187-197.
- Zhao, Z.-S., Leung, T., Manser, E., and Lim, L. 1995. Pheromone signalling in *Saccharomyces cerevisiae* requires the small GTP-binding protein Cdc42p and its activator CDC24. Mol. Cell. Biol. 15: 5246-5257.
- Simon, M.-N., De Virgilio, C., Souza, B., Pringle, J.R., Abo, A., and Reed, S.I. 1995. Role for the Rho-family GTPase Cdc42 in yeast matingpheromone signal pathway. Nature 376: 702-705.
- Bardwell, L., Cook, J.G., Chang, E.C., Cairns, B.R., and Thorner, J. 1996. Signaling in the yeast pheromone response pathway: specific and highaffinity interactions of the mitogen-activated protein (MAP) kinases Kss1 and Fus3 with the upstream MAP kinase kinase Ste7. Mol. Cell. Biol. 16: 3637-3650.

SOURCE

Cdc24 (yN-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of Cdc24 of *Saccharomyces cerevisiae* origin.

PRODUCT

Each vial contains 200 μg lgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-6794 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Cdc24 (yN-18) is recommended for detection of Cdc24 of *Saccharomyces cerevisiae* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2033 and Western Blotting Luminol Reagent: sc-2048.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.