

PLSCR3 (T-13): sc-68123

BACKGROUND

The calcium-dependent mitochondrial membrane protein PLSCR3 (phospholipid scramblase 3), also known as PLS3, is a member of the phospholipid scramblase (PLS) family. The PLS family consists of membrane-bound enzymes that participate in the bi-directional movement of phospholipids. PLSCR3 is expressed in a wide variety of tissues but is not found in brain, liver or testis. It is involved in the regulation of mitochondrial respiratory function, morphology and apoptotic responses. More specifically, PLSCR3 is responsible for mediating the transport of cardiolipin from the inner mitochondrial membrane to the outer mitochondrial membrane. Cardiolipin is a major polyglycerophospholipid that plays a role in the regulation of mitochondrial enzymes involved in the generation of ATP. PLSCR3 activity is activated/phosphorylated by PKC δ and PLSCR3 functions as a downstream effector in PKC δ -induced apoptosis.

REFERENCES

1. Wiedmer, T., Zhou, Q., Kwok, D.Y. and Sims, P.J. 2000. Identification of three new members of the phospholipid scramblase gene family. *Biochim. Biophys. Acta* 1467: 244-253.
2. Liu, J., Dai, Q., Chen, J., Durrant, D., Freeman, A., Liu, T., Grossman, D. and Lee, R.M. 2003. Phospholipid scramblase 3 controls mitochondrial structure, function, and apoptotic response. *Mol. Cancer Res.* 1: 892-902.
3. Liu, J., Chen, J., Dai, Q. and Lee, R.M. 2003. Phospholipid scramblase 3 is the mitochondrial target of protein kinase C δ -induced apoptosis. *Cancer Res.* 63: 1153-1156.
4. Wiedmer, T., Zhao, J., Li, L., Zhou, Q., Hevener, A., Olefsky, J.M., Curtiss, L.K. and Sims, P.J. 2004. Adiposity, dyslipidemia, and Insulin resistance in mice with targeted deletion of phospholipid scramblase 3 (PLSCR3). *Proc. Natl. Acad. Sci. USA* 101: 13296-13301.
5. He, Y., Liu, J., Durrant, D., Yang, H.S., Sweatman, T., Lothstein, L. and Lee, R.M. 2005. N-benzyladriamycin-14-valerate (AD198) induces apoptosis through protein kinase C- δ -induced phosphorylation of phospholipid scramblase 3. *Cancer Res.* 65: 10016-10023.
6. Phillippe, M., Bradley, D.F., Ji, H., Oppenheimer, K.H. and Chien, E.K. 2006. Phospholipid scramblase isoform expression in pregnant rat uterus. *J. Soc. Gynecol. Investig.* 13: 497-501.
7. Van, Q., Liu, J., Lu, B., Feingold, K.R., Shi, Y., Lee, R.M. and Hatch, G.M. 2007. Phospholipid scramblase-3 regulates cardiolipin *de novo* biosynthesis and its resynthesis in growing HeLa cells. *Biochem. J.* 401: 103-109.
8. Mutch, D.M., O'Maille, G., Wikoff, W.R., Wiedmer, T., Sims, P.J. and Siuzdak, G. 2007. Mobilization of pro-inflammatory lipids in obese PLSCR3-deficient mice. *Genome Biol.* 8: R38.

CHROMOSOMAL LOCATION

Genetic locus: PLSCR3 (human) mapping to 17p13.1; Plscr3 (mouse) mapping to 11 B3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

PLSCR3 (T-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of PLSCR3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68123 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PLSCR3 (T-13) is recommended for detection of PLSCR3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PLSCR3 (T-13) is also recommended for detection of PLSCR3 in additional species, including equine and bovine.

Suitable for use as control antibody for PLSCR3 siRNA (h): sc-62828, PLSCR3 siRNA (m): sc-62829, PLSCR3 shRNA Plasmid (h): sc-62828-SH, PLSCR3 shRNA Plasmid (m): sc-62829-SH, PLSCR3 shRNA (h) Lentiviral Particles: sc-62828-V and PLSCR3 shRNA (m) Lentiviral Particles: sc-62829-V.

Molecular Weight of PLSCR3: 32 kDa.

Positive Controls: 3T3-L1 cell lysate: sc-2243, Sol8 cell lysate: sc-2249 or mouse lung extract: sc-2390.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.