

PLSCR4 (N-12): sc-68126

BACKGROUND

The calcium-dependent plasma membrane protein PLSCR4 (phospholipid scramblase 4), also known as PLS4, TRA1 or cell growth inhibiting protein 43, is a 329 amino acid member of the phospholipid scramblase (PLS) family. The PLS family consists of membrane-bound enzymes that participate in the bidirectional movement of phospholipids. PLSCR4 is expressed in a wide variety of tissues including brain (unlike other PLS family members), but it is not found in peripheral blood lymphocytes. It contains an N-terminal proline-rich region with several PxxP motifs and a calcium-binding domain at its C-terminus. PLSCR4 shares 46% identity with PLSCR1.

REFERENCES

1. Wiedmer, T., Zhou, Q., Kwok, D.Y. and Sims, P.J. 2000. Identification of three new members of the phospholipid scramblase gene family. *Biochim. Biophys. Acta* 1467: 244-253.
2. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607612. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Wiedmer, T., Zhao, J., Li, L., Zhou, Q., Hevener, A., Olefsky, J.M., Curtiss, L.K. and Sims, P.J. 2004. Adiposity, dyslipidemia, and Insulin resistance in mice with targeted deletion of phospholipid scramblase 3 (PLSCR3). *Proc. Natl. Acad. Sci. USA* 101: 13296-13301.
4. Chen, M.H., Ben-Efraim, I., Mitrousis, G., Walker-Kopp, N., Sims, P.J. and Cingolani, G. 2005. Phospholipid scramblase 1 contains a nonclassical nuclear localization signal with unique binding site in Importin α . *J. Biol. Chem.* 280: 10599-10606.
5. Phillippe, M., Bradley, D.F., Ji, H., Oppenheimer, K.H. and Chien, E.K. 2006. Phospholipid scramblase isoform expression in pregnant rat uterus. *J. Soc. Gynecol. Investig.* 13: 497-501.
6. Sahu, S.K., Gummadi, S.N., Manoj, N. and Aradhyam, G.K. 2007. Phospholipid scramblases: an overview. *Arch. Biochem. Biophys.* 462: 103-114.
7. Lu, B., Sims, P.J., Wiedmer, T., Moser, A.H., Shigenaga, J.K., Grunfeld, C. and Feingold, K.R. 2007. Expression of the phospholipid scramblase (PLSCR) gene family during the acute phase response. *Biochim. Biophys. Acta* 1771: 1177-1185.
8. Kim, S., Choi, K.H., Baykiz, A.F. and Gershenfeld, H.K. 2007. Suicide candidate genes associated with bipolar disorder and schizophrenia: an exploratory gene expression profiling analysis of post-mortem prefrontal cortex. *BMC Genomics* 8: 413-413.

CHROMOSOMAL LOCATION

Genetic locus: PLSCR4 (human) mapping to 3q24; Plscr4 (mouse) mapping to 9 E3.3.

SOURCE

PLSCR4 (N-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of PLSCR4 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68126 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PLSCR4 (N-12) is recommended for detection of phospholipid scramblase 4 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PLSCR4 (N-12) is also recommended for detection of phospholipid scramblase 4 in additional species, including equine and canine.

Suitable for use as control antibody for PLSCR4 siRNA (h): sc-62830, PLSCR4 siRNA (m): sc-62831, PLSCR4 shRNA Plasmid (h): sc-62830-SH, PLSCR4 shRNA Plasmid (m): sc-62831-SH, PLSCR4 shRNA (h) Lentiviral Particles: sc-62830-V and PLSCR4 shRNA (m) Lentiviral Particles: sc-62831-V.

Molecular Weight of PLSCR4: 37 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.