

frizzled-2 (S-13): sc-68330

BACKGROUND

The frizzled gene, originally identified in *Drosophila melanogaster*, is involved in the development of tissue polarity. The frizzled proteins contain seven transmembrane domains, a cysteine-rich domain in the extracellular region and a carboxy-terminal Ser/Thr-XXX-Val motif. The proteins in this family function as receptors for Wnt and are generally coupled to G proteins. Frizzleds are members of the G protein-coupled receptor superfamily. Frizzled-2 is expressed in the fetal kidney and lung and in the adult ovary and colon. Frizzled-2 mediates the Wnt/cGMP/Ca²⁺ pathway. It binds Wnt proteins and signals by activating the release of stored calcium. frizzled-2 expression is regulated by Angiotensin II. Activated frizzled-2 suppresses the activity of protein kinase G, and activates NFAT-dependent transcription, the phosphatidylinositol pathway and calcium sensitive enzymes.

REFERENCES

- Zhao, Z., et al. 1995. A human homologue of the *Drosophila* polarity gene frizzled has been identified and mapped to 17q21.1. *Genomics* 27: 370-373.
- Sagara, N., et al. 1998. Molecular cloning, differential expression, and chromosomal localization of human frizzled-1, frizzled-2, and frizzled-7. *Biochem. Biophys. Res. Commun.* 252: 117-122.
- Malbon, C.C., et al. 2004. Frizzleds: new members of the superfamily of G protein-coupled receptors. *Front. Biosci.* 9: 1048-1058.
- Castoldi, G., et al. 2005. Angiotensin II modulates frizzled-2 receptor expression in rat vascular smooth muscle cells. *Clin. Sci.* 108: 523-530.
- Rodriguez, J., et al. 2005. SFRP1 regulates the growth of retinal ganglion cell axons through the Fz2 receptor. *Nat. Neurosci.* 8: 1301-1309.
- Jiang, F., et al. 2006. Gene expression profile of quiescent and activated rat hepatic stellate cells implicates Wnt signaling pathway in activation. *J. Hepatol.* 45: 401-409.
- Ma, L., et al. 2006. Suppression of cyclic GMP-dependent protein kinase is essential to the Wnt/cGMP/Ca²⁺ pathway. *J. Biol. Chem.* 281: 30990-31001.
- Harris, K.E., et al. 2007. Different Wnt signals act through the frizzled and RYK receptors during *Drosophila* salivary gland migration. *Development* 134: 2017-2025.
- Hurson, C.J., et al. 2007. Gene expression analysis in human osteoblasts exposed to dexamethasone identifies altered developmental pathways as putative drivers of osteoporosis. *BMC Musculoskelet. Disord.* 8: 12.

CHROMOSOMAL LOCATION

Genetic locus: FZD2 (human) mapping to 17q21.31; Fzd2 (mouse) mapping to 11 E1.

SOURCE

frizzled-2 (S-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an extracellular domain of frizzled-2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68330 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

frizzled-2 (S-13) is recommended for detection of frizzled-2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

frizzled-2 (S-13) is also recommended for detection of frizzled-2 in additional species, including canine and porcine.

Suitable for use as control antibody for frizzled-2 siRNA (h): sc-39979, frizzled-2 siRNA (m): sc-39980, frizzled-2 shRNA Plasmid (h): sc-39979-SH, frizzled-2 shRNA Plasmid (m): sc-39980-SH, frizzled-2 shRNA (h) Lentiviral Particles: sc-39979-V and frizzled-2 shRNA (m) Lentiviral Particles: sc-39980-V.

Molecular Weight of frizzled-2: 65 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **Bfrizzled-2 (K19R): sc-74019**, our highly recommended monoclonal alternative to frizzled-2 (S-13).