

Inhibin β -E (E-11): sc-68566

BACKGROUND

Inhibin β -E, also known as INHBE, is a 350 amino acid secreted protein that belongs to the TGF β superfamily and exists as either a homodimer or a heterodimer via association with other Inhibin subunits. Functioning to inhibit the secretion of follitropin by the pituitary gland, Inhibin β -E is involved in a variety of events throughout the body, including germ cell development, nerve cell survival, insulin secretion and bone growth. The gene encoding Inhibin β -E maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

- Allen, T.L., et al. 1996. Cytogenetic and molecular analysis in trisomy 12p. *Am. J. Med. Genet.* 63: 250-256.
- Yang, W., et al. 1998. Low basal transcripts of the COL2A1 collagen gene from lymphoblasts show alternative splicing of exon 12 in the Kniest form of spondyloepiphyseal dysplasia. *Hum. Mutat.* 1: S1-S2.
- Trowsdale, J., et al. 2001. The genomic context of natural killer receptor extended gene families. *Immunol. Rev.* 181: 20-38.
- Hashimoto, O., et al. 2002. cDNA cloning and expression of human activin β E subunit. *Mol. Cell. Endocrinol.* 194: 117-122.
- Zumkeller, W., et al. 2004. Genotype/phenotype analysis in a patient with pure and complete trisomy 12p. *Am. J. Med. Genet. A* 129: 261-264.
- Kelley, J., et al. 2005. Comparative genomics of natural killer cell receptor gene clusters. *PLoS Genet.* 1: e27.
- Segel, R., et al. 2006. The natural history of trisomy 12p. *Am. J. Med. Genet. A* 140: 695-703.
- Hashimoto, O., et al. 2006. Impaired growth of pancreatic exocrine cells in transgenic mice expressing human activin β E subunit. *Biochem. Biophys. Res. Commun.* 341: 416-424.
- Online Mendelian Inheritance in Man, OMIM™. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 612031. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: INHBE (human) mapping to 12q13.3; Inhbe (mouse) mapping to 10 D3.

SOURCE

Inhibin β -E (E-11) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Inhibin β -E of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68566 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Inhibin β -E (E-11) is recommended for detection of Inhibin β -E of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Inhibin β -E (E-11) is also recommended for detection of Inhibin β -E in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Inhibin β -E siRNA (h): sc-75338, Inhibin β -E siRNA (m): sc-75339, Inhibin β -E shRNA Plasmid (h): sc-75338-SH, Inhibin β -E shRNA Plasmid (m): sc-75339-SH, Inhibin β -E shRNA (h) Lentiviral Particles: sc-75338-V and Inhibin β -E shRNA (m) Lentiviral Particles: sc-75339-V.

Molecular Weight of Inhibin β -E: 39 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.