

PKDREJ (Y-15): sc-68725

BACKGROUND

PKDREJ (polycystic kidney disease and receptor for egg jelly-related protein) is a 2,253 amino acid multi-pass membrane protein that belongs to the polycystin family and is encoded by an intronless gene. Expressed exclusively in testis, PKDREJ is thought to play a key role in fertilization, possibly by generating a calcium transport channel that is directly associated with the acrosome reaction of the sperm. Human PKDREJ, which can exist as homomultimers or heteromultimers, shares 64% identity with its mouse homolog, suggesting a conserved function between species. PKDREJ contains one receptor for egg jelly (REJ) domain, eleven transmembrane domains, one G-protein-coupled receptor proteolytic site (GPS) domain and one lipoygenase domain. Several isoforms of PKDREJ exist due to alternative splicing events.

REFERENCES

- Moy, G.W., Mendoza, L.M., Schulz, J.R., Swanson, W.J., Glabe, C.G. and Vacquier, V.D. 1996. The sea urchin sperm receptor for egg jelly is a modular protein with extensive homology to the human polycystic kidney disease protein, PKD1. *J. Cell Biol.* 133: 809-817.
- Hughes, J., Ward, C.J., Aspinwall, R., Butler, R. and Harris, P.C. 1999. Identification of a human homologue of the sea urchin receptor for egg jelly: a polycystic kidney disease-like protein. *Hum. Mol. Genet.* 8: 543-549.
- Mengerink, K.J., Moy, G.W. and Vacquier, V.D. 2000. suREJ proteins: new signalling molecules in sea urchin spermatozoa. *Zygote* 8: 28-30.
- Qian, F., Boletta, A., Bhunia, A.K., Xu, H., Liu, L., Ahrabi, A.K., Watnick, T.J., Zhou, F. and Germino, G.G. 2002. Cleavage of polycystin-1 requires the receptor for egg jelly domain and is disrupted by human autosomal-dominant polycystic kidney disease 1-associated mutations. *Proc. Natl. Acad. Sci. USA* 99: 16981-16986.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604670. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Galindo, B.E., Moy, G.W. and Vacquier, V.D. 2004. A third sea urchin sperm receptor for egg jelly module protein, suREJ2, concentrates in the plasma membrane over the sperm mitochondrion. *Dev. Growth Differ.* 46: 53-60.
- Sutton, K.A., Jungnickel, M.K., Ward, C.J., Harris, P.C. and Florman, H.M. 2006. Functional characterization of PKDREJ, a male germ cell-restricted polycystin. *J. Cell. Physiol.* 209: 493-500.
- Butscheid, Y., Chubanov, V., Steger, K., Meyer, D., Dietrich, A. and Gudermann, T. 2006. Polycystic kidney disease and receptor for egg jelly is a plasma membrane protein of mouse sperm head. *Mol. Reprod. Dev.* 73: 350-360.
- Hamm, D., Mautz, B.S., Wolfner, M.F., Aquadro, C.F. and Swanson, W.J. 2007. Evidence of amino acid diversity-enhancing selection within humans and among primates at the candidate sperm-receptor gene PKDREJ. *Am. J. Hum. Genet.* 81: 44-52.

CHROMOSOMAL LOCATION

Genetic locus: PKDREJ (human) mapping to 22q13.31; Pkdrej (mouse) mapping to 15 E2.

SOURCE

PKDREJ (Y-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of PKDREJ of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68725 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PKDREJ (Y-15) is recommended for detection of PKDREJ of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for PKDREJ siRNA (h): sc-76157, PKDREJ siRNA (m): sc-76158, PKDREJ shRNA Plasmid (h): sc-76157-SH, PKDREJ shRNA Plasmid (m): sc-76158-SH, PKDREJ shRNA (h) Lentiviral Particles: sc-76157-V and PKDREJ shRNA (m) Lentiviral Particles: sc-76158-V.

Molecular Weight of PKDREJ: 241 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.