

SMPDL3B (M-76): sc-68844

BACKGROUND

SMPDL3B (sphingomyelin phosphodiesterase, acid-like 3B), also known as ASMLPD or ASML3B (acid sphingomyelinase-like phosphodiesterase 3B), is a 455 amino acid secreted protein belonging to the acid sphingomyelinase family. It is expressed in granulosa cells of the ovarian follicle and is a homolog of ASM (acid sphingomyelinase). ASM is a lysosomal protein that hydrolyzes sphingomyelin to ceramide and phosphocholine playing an important role in apoptosis of germ cell lines. Deficiency of ASM is associated with type A and type B Niemann-Pick disease. Type A is a fatal neurodegenerative disorder seen in infancy and resulting in death by age three, whereas type B is a non-neuropathic disease with a later onset.

REFERENCES

- Quintern, L.E., Weitz, G., Nehrkorn, H., Tager, J.M., Schram, A.W. and Sandhoff, K. 1987. Acid sphingomyelinase from human urine: purification and characterization. *Biochim. Biophys. Acta* 922: 323-336.
- Schuchman, E.H., Suchi, M., Takahashi, T., Sandhoff, K. and Desnick, R.J. 1991. Human acid sphingomyelinase. Isolation, nucleotide sequence and expression of the full-length and alternatively spliced cDNAs. *J. Biol. Chem.* 266: 8531-8539.
- Levrán, O., Desnick, R.J. and Schuchman, E.H. 1991. Niemann-Pick disease: a frequent missense mutation in the acid sphingomyelinase gene of Ashkenazi Jewish type A and B patients. *Proc. Natl. Acad. Sci. USA* 88: 3748-3752.
- Takahashi, T., Suchi, M., Desnick, R.J., Takada, G. and Schuchman, E.H. 1992. Identification and expression of five mutations in the human acid sphingomyelinase gene causing types A and B Niemann-Pick disease. Molecular evidence for genetic heterogeneity in the neuronopathic and non-neuronopathic forms. *J. Biol. Chem.* 267: 12552-12558.
- Langmann, T., Buechler, C., Ries, S., Schaeffler, A., Aslanidis, C., Schuierer, M., Weiler, M., Sandhoff, K., de Jong, P.J. and Schmitz, G. 1999. Transcription factors Sp1 and AP-2 mediate induction of acid sphingomyelinase during monocytic differentiation. *J. Lipid Res.* 40: 870-880.
- Perosa, F., Favoino, E., Caragnano, M.A. and Dammacco, F. 2006. Generation of biologically active linear and cyclic peptides has revealed a unique fine specificity of rituximab and its possible cross-reactivity with acid sphingomyelinase-like phosphodiesterase 3b precursor. *Blood* 107: 1070-1077.
- Rizzolio, F., Sala, C., Alboresi, S., Bione, S., Gilli, S., Goegan, M., Pramparo, T., Zuffardi, O. and Toniolo, D. 2007. Epigenetic control of the critical region for premature ovarian failure on autosomal genes translocated to the X chromosome: a hypothesis. *Hum. Genet.* 121: 441-450.

CHROMOSOMAL LOCATION

Genetic locus: SMPDL3B (human) mapping to 1p35.3; Smpdl3b (mouse) mapping to 4 D2.3.

SOURCE

SMPDL3B (M-76) is a rabbit polyclonal antibody raised against amino acids 196-271 mapping within an internal region of SMPDL3B of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

SMPDL3B (M-76) is recommended for detection of SMPDL3B of mouse, rat and, to a lesser extent, human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1–2 µg per 100–500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Suitable for use as control antibody for SMPDL3B siRNA (h): sc-76525, SMPDL3B siRNA (m): sc-76526, SMPDL3B shRNA Plasmid (h): sc-76525-SH, SMPDL3B shRNA Plasmid (m): sc-76526-SH, SMPDL3B shRNA (h) Lentiviral Particles: sc-76525-V and SMPDL3B shRNA (m) Lentiviral Particles: sc-76526-V.

Molecular Weight of SMPDL3B: 51 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



SMPDL3B (M-76): sc-68844. Western blot analysis of SMPDL3B expression in non-transfected: sc-117752 (A) and mouse SMPDL3B transfected: sc-126022 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.