

Bag-2 (E-20): sc-68949

BACKGROUND

Bag-2 (Bcl-2-associated athanogene 2), also known as Bag family molecular chaperone regulator 2, is a member of the Bag family of proteins and contains the most diverged of the characteristic C-terminal Bag domain. Via their Bag domain, Bag proteins bind with high affinity to the HSC 70/HSP 70 ATPase domain, regulating chaperone activity and apoptosis. Bag-2 is an evolutionarily conserved cytoplasmic protein with putative N-terminal phosphorylation sites and specifically functions as an HSC 70 co-chaperone. Bag-2 is a major component of the HSC 70/CHIP chaperone-dependent ubiquitin ligase complex and acts to disrupt CHIP-mediated ubiquitylation. In this complex, Bag-2 directly interacts with the ATPase domain of HSC 70 as well as the U-box domain of CHIP and inhibits ubiquitylation by interfering with the association between CHIP and its ubiquitin conjugating enzyme.

REFERENCES

1. Takayama, S., Xie, Z. and Reed, J.C. 1999. An evolutionarily conserved family of HSP70/HSC70 molecular chaperone regulators. *J. Biol. Chem.* 274: 781-786.
2. Ueda, K., Kosako, H., Fukui, Y. and Hattori, S. 2004. Proteomic identification of Bcl2-associated athanogene 2 as a novel MAPK-activated protein kinase 2 substrate. *J. Biol. Chem.* 279: 41815-41821.
3. Arndt, V., Daniel, C., Nastainczyk, W., Alberti, S. and Höhfeld, J. 2005. Bag-2 acts as an inhibitor of the chaperone-associated ubiquitin ligase CHIP. *Mol. Biol. Cell* 16: 5891-5900.
4. Dai, Q., Qian, S.B., Li, H.H., McDonough, H., Borchers, C., Huang, D., Takayama, S., Younger, J.M., Ren, H.Y., Cyr, D.M. and Patterson, C. 2005. Regulation of the cytoplasmic quality control protein degradation pathway by Bag-2. *J. Biol. Chem.* 280: 38673-38681.
5. Götz, R., Wiese, S., Takayama, S., Camarero, G.C., Rossoll, W., Schweizer, U., Troppmair, J., Jablonka, S., Holtmann, B., Reed, J.C., Rapp, U.R. and Sendtner, M. 2005. Bag-1 is essential for differentiation and survival of hematopoietic and neuronal cells. *Nat. Neurosci.* 8: 1169-1178.
6. Wada, S., Hamada, M. and Satoh, N. 2006. A genomewide analysis of genes for the heat shock protein 70 chaperone system in the ascidian *Ciona intestinalis*. *Cell Stress Chaperones* 11: 23-33.
7. Yi, C., Li, S., Wang, J., Wei, N. and Deng, X.W. 2006. Affinity purification reveals the association of WD40 protein constitutive photomorphogenic 1 with the hetero-oligomeric TCP-1 chaperonin complex in mammalian cells. *Int. J. Biochem. Cell Biol.* 38: 1076-1083.
8. Arndt, V., Rogon, C. and Höhfeld, J. 2007. To be, or not to be—molecular chaperones in protein degradation. *Cell. Mol. Life Sci.* 64: 2525-2541.

CHROMOSOMAL LOCATION

Genetic locus: BAG2 (human) mapping to 6p11.2; Bag2 (mouse) mapping to 1 B.

SOURCE

Bag-2 (E-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Bag-2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-68949 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Bag-2 (E-20) is recommended for detection of Bag-2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Bag-2 (E-20) is also recommended for detection of Bag-2 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for Bag-2 siRNA (h): sc-72600, Bag-2 siRNA (m): sc-72601, Bag-2 shRNA Plasmid (h): sc-72600-SH, Bag-2 shRNA Plasmid (m): sc-72601-SH, Bag-2 shRNA (h) Lentiviral Particles: sc-72600-V and Bag-2 shRNA (m) Lentiviral Particles: sc-72601-V.

Molecular Weight of Bag-2: 26 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200 or Jurkat whole cell lysate: sc-2204.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **Bag-2 (C-6): sc-390107** or **Bag-2 (A-7): sc-390262**, our highly recommended monoclonal alternatives to Bag-2 (E-20).