

Fra-1 (C-14): sc-69271

BACKGROUND

The v-Fos oncogene was initially detected in two independent murine osteosarcoma virus isolates and an avian nephroblastoma virus. Members of the c-Fos gene family, including c-Fos, Fos B, Fra-1 and Fra-2, encode nuclear phosphoproteins that are rapidly and transiently induced by a variety of agents and function as transcriptional regulators for several genes. In contrast to c-Jun proteins which form homo- and heterodimers which bind to specific DNA response elements, c-Fos proteins are only active as heterodimers with members of the Jun gene family. In addition, selected ATF/CREB family members can form leucine zipper dimers with Fos and Jun. Different dimers exhibit differential specificity and affinity for AP-1 and CRE sites.

REFERENCES

1. Finkel, M.P., Biskis, B.O. and Jenkins, P.B. 1966. Virus induction of osteosarcomas in mice. *Science* 151: 698-701.
2. Sambucetti, L.C. and Curran, T. 1986. The Fos protein complex is associated with DNA in isolated nuclei and binds to DNA cellulose. *Science* 234: 1417-1419.
3. Nishizawa, M., Goto, N. and Kawai, S. 1987. An avian transforming retrovirus isolated from a nephroblastoma that carries the Fos gene as the oncogene. *J. Virol.* 61: 3733-3740.
4. Bohmann, D., Bos, T.J., Admon, A., Nishimura, T., Vogt, P.K. and Tjian, R. 1987. Human proto-oncogene c-Jun encodes a DNA binding protein with structural and functional properties of transcription factor AP-1. *Science* 238: 1386-1392.
5. Cohen, D.R., Ferreira, P.C., Gentz, R., Franza, B.R., Jr. and Curran, T. 1989. The product of a Fos-related gene, Fra-1, binds cooperatively to the AP-1 site with Jun: transcription factor AP-1 is comprised of multiple protein complexes. *Genes Dev.* 3: 173-184.
6. Nishina, H., Sato, H., Suzuki, T., Sato, M. and Iba, H. 1990. Isolation and characterization of Fra-2, an additional member of the Fos gene family. *Proc. Natl. Acad. Sci. USA* 87: 3619-3623.
7. Boise, L.H., Petryniak, B., Mao, X., June, C.H., Wang, C.Y., Lindsten, T., Bravo, R., Kovary, K., Leiden, J.M. and Thompson, C.B. 1993. The NFAT-1 DNA binding complex in activated T cells contains Fra-1 and Jun B. *Mol. Cell. Biol.* 13: 1911-1919.
8. Casalino, L., Bakiri, L., Talotta, F., Weitzman, J.B., Fusco, A., Yaniv, M. and Verde, P. 2007. Fra-1 promotes growth and survival in Ras-transformed thyroid cells by controlling cyclin A transcription. *EMBO J.* 26: 1878-1890.
9. Chiappetta, G., Ferraro, A., Botti, G., Monaco, M., Pasquinelli, R., Vuttariello, E., Arnaldi, L., Di Bonito, M., D'Aiuto, G., Pierantoni, G.M. and Fusco, A. 2007. Fra-1 protein overexpression is a feature of hyperplastic and neoplastic breast disorders. *BMC Cancer* 7: 17.

CHROMOSOMAL LOCATION

Genetic locus: Fos1 (mouse) mapping to 19 A.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

Fra-1 (C-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Fra-1 of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-69271 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-69271 X, 200 µg/0.1 ml.

APPLICATIONS

Fra-1 (C-14) is recommended for detection of Fra-1 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Fra-1 siRNA (m): sc-35406, Fra-1 shRNA Plasmid (m): sc-35406-SH and Fra-1 shRNA (m) Lentiviral Particles: sc-35406-V.

Fra-1 (C-14) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of Fra-1: 40 kDa.

Positive Controls: NIH/3T3 nuclear extract: sc-2138, 3611-RF nuclear extract: sc-2143 or 3611-RF + PMA nuclear extract: sc-2144.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.



Try **Fra-1 (D-3): sc-376148** or **Fra-1 (C-12): sc-28310**, our highly recommended monoclonal alternatives to Fra-1 (C-14). Also, for AC, HRP, FITC, PE, Alexa Fluor[®] 488 and Alexa Fluor[®] 647 conjugates, see **Fra-1 (D-3): sc-376148**.