HIV-1 Integrase (IN-2): sc-69721



The Power to Question

BACKGROUND

The human immunodeficiency virus type 1 (HIV-1) uses an enzyme HIV-1 Integrase to carry out the integration of its viral DNA into the host chromosome, thereby tricking the host cell machinery into making viral proteins. HIV-1 Integrase is a viral protein that catalyzes the first two steps of the three-step DNA integration process and belongs to the superfamily of polynucleotidyl transferases. The integration of HIV-1 DNA into the host chromosome is achieved by the integrase protein performing a series of DNA cleaving and ligation reactions. The C-terminal domain of HIV-1 Integrase forms a homodimer, consisting of two structures that are similar to Src-homology 3 (SH3) domains, while the N-terminal domain of HIV-1 Integrase contains a pair of His and Cys residues (the HHCC motif) that are conserved among retroviral integrases. Zinc binds to the HHCC motif and supports the folded state of the N-terminus.

REFERENCES

- Engelman, A., et al. 1993. Identification of discrete functional domains of HIV-1 Integrase and their organization within an active multimeric complex. EMBO J. 12: 3269-3275.
- Dyda, F., et al. 1994. Crystal structure of the catalytic domain of HIV-1 Integrase: similarity to other polynucleotidyl transferases. Science 266: 1981-1986.
- 3. Eijkelenboom, A.P., et al. 1995. The DNA-binding domain of HIV-1 Integrase has an SH3-like fold. Nat. Struct. Biol. 2: 807-810.

SOURCE

 $\mbox{HIV-1}$ Integrase (IN-2) is a mouse monoclonal antibody raised against $\mbox{HIV-1}$ Integrase.

PRODUCT

Each vial contains 200 $\mu g \ lgG_1$ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

HIV-1 Integrase (IN-2) is available conjugated to agarose (sc-69721 AC), 500 $\mu g/0.25$ ml agarose in 1 ml, for IP; to HRP (sc-69721 HRP), 200 $\mu g/ml$, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-69721 PE), fluorescein (sc-69721 FITC), Alexa Fluor® 488 (sc-69721 AF488), Alexa Fluor® 546 (sc-69721 AF546), Alexa Fluor® 594 (sc-69721 AF594) or Alexa Fluor® 647 (sc-69721 AF647), 200 $\mu g/ml$, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-69721 AF680) or Alexa Fluor® 790 (sc-69721 AF790), 200 $\mu g/ml$, for Near-Infrared (NIR) WB, IF and FCM.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA

APPLICATIONS


HIV-1 Integrase (IN-2) is recommended for detection of Pol Integrase of HIV-1 by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)] and immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

Molecular Weight of HIV-1 Integrase: 32 kDa.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

DATA



HIV-1 Integrase (IN-2): sc-69721. Western blot analysis of HIV-1 Integrase expression in semi-purified HIV-1

SELECT PRODUCT CITATIONS

- Thierry, S., et al. 2015. Integrase inhibitor reversal dynamics indicate unintegrated HIV-1 DNA initiate de novo integration. Retrovirology 12: 24.
- Zhang, X., et al. 2016. Identification of interaction between HIV-1 glycoprotein 41 and integrase. Virol. Sin. 31: 415-424.
- Martinez, Z.S., et al. 2016. Fullerene derivatives strongly inhibit HIV-1 replication by affecting virus maturation without impairing protease activity. Antimicrob. Agents Chemother. 60: 5731-5741.
- Kim, H., et al. 2018. HBV polymerase-derived peptide exerts an anti-HIV-1 effect by inhibiting the acetylation of viral integrase. Biochem. Biophys. Res. Commun. 501: 541-546.
- 5. Ali, H., et al. 2019. Cellular TRIM33 restrains HIV-1 infection by targeting viral integrase for proteasomal degradation. Nat. Commun. 10: 926.
- 6. Yadav, P., et al. 2019. Interaction of HIV-1 integrase with polypyrimidine tract binding protein and associated splicing factor (PSF) and its impact on HIV-1 replication. Retrovirology 16: 12.
- Allouch, A., et al. 2020. SUGT1 controls susceptibility to HIV-1 infection by stabilizing microtubule plus-ends. Cell Death Differ. 27: 3243-3257.
- 8. Mete, B., et al. 2022. Human immunodeficiency virus type 1 impairs sumoylation. Life Sci. Alliance 5: e202101103.
- 9. Liang, G., et al. 2022. CTNNBL1 restricts HIV-1 replication by suppressing viral DNA integration into the cell genome. Cell Rep. 38: 110533.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

Santa Cruz Biotechnology, Inc. 1.800.457.3801 831.457.3801 fax 831.457.3801 Europe +00800 4573 8000 49 6221 4503 0 www.scbt.com