

Hrk (M-20): sc-6973

BACKGROUND

Members of the Bcl-2 family of proteins interact to regulate programmed cell death, or apoptosis. Various homodimers and heterodimers formed by proteins in this family can either promote or inhibit apoptosis. Bcl-2 blocks cell death following a variety of stimuli and confers a death-sparing effect on certain hematopoietic cell lines following growth factor withdrawal. Additional apoptotic inhibitors in this family include Bcl-x, Bcl-w, Mcl-1, Bag-1 and A1. Pro-apoptotic members of this family include Bax, Bad, Bak, NBK (Bik), BID and Hrk. Hrk (for harakiri), designated DP5 or neuronal death protein in mouse and rat, contains a BH3 domain with high homology to other Bcl-2 family members but lacks the conserved BH1 and BH2 domains. Physical interaction of Hrk with Bcl-2 or Bcl-x_L inhibits the apoptotic activity of Hrk.

REFERENCES

- Vaux, D.L., et al. 1988. Bcl-2 promotes the survival of hemopoietic cells and cooperates with c-Myc to immortalize pre-B cells. *Nature* 335: 440-442.
- Nunez, G., et al. 1990. Deregulated Bcl-2 gene expression selectively prolongs survival of growth factor-deprived hemopoietic cell lines. *J. Immunol.* 144: 3602-3610.
- Oltvai, Z.N., et al. 1993. Bcl-2 heterodimerizes *in vivo* with a conserved homolog, Bax, that accelerates programmed cell death. *Cell* 74: 609-619.
- Sato, T., et al. 1994. Interactions among members of the Bcl-2 protein family analyzed with a yeast two-hybrid system. *Proc. Natl. Acad. Sci. USA* 91: 9238-9242.
- Oltvai, Z.N., et al. 1994. Checkpoints of dueling dimers foil death wishes. *Cell* 79: 189-192.
- Yang, E., et al. 1996. Molecular thanatopsis: a discourse on the Bcl-2 family and cell death. *Blood* 88: 386-401.
- Wang, K., et al. 1996. BID: a novel BH3 domain-only death agonist. *Genes Dev.* 10: 2859-2869.
- Nagata, S. 1997. Apoptosis by death factor. *Cell* 88: 355-365.
- Inohara, N., et al. 1997. Harakiri, a novel regulator of cell death, encodes a protein that activates apoptosis and interacts selectively with survival-promoting proteins Bcl-2 and Bcl-x_L. *EMBO J.* 16: 1686-1694.

CHROMOSOMAL LOCATION

Genetic locus: Bid3 (mouse) mapping to 5 F.

SOURCE

Hrk (M-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Hrk of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-6973 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Hrk (M-20) is recommended for detection of Hrk of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Hrk siRNA (m): sc-37304, Hrk shRNA Plasmid (m): sc-37304-SH and Hrk shRNA (m) Lentiviral Particles: sc-37304-V.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

SELECT PRODUCT CITATIONS

- Korhonen, L., et al. 2003. Increase in Bcl-2 phosphorylation and reduced levels of BH3-only Bcl-2 family proteins in kainic acid-mediated neuronal death in the rat brain. *Eur. J. Neurosci.* 18: 1121-1134.
- Cartron, P.F., et al. 2012. Prognostic impact of the expression/phosphorylation of the BH3-only proteins of the Bcl-2 family in glioblastoma multiforme. *Cell Death Dis.* 3: e421.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.