

Che-1 (C-19): sc-70032

BACKGROUND

Che-1 (protein AATF, apoptosis-antagonizing transcription factor) is a widely expressed nuclear protein that belongs to the AATF family. Hyperphosphorylated during the G₁/S phase transition, Che-1 may function as a general inhibitor of the histone deacetylase HDAC1. Che-1 binding to the pocket region of Rb may displace HDAC1 from Rb/E2F complexes, leading to activation of E2F target genes and cell cycle progression. Displacement of HDAC1 from Sp1 bound to the p21 promoter leads to increased expression of Che-1. It also antagonizes PAR4 (prostate apoptosis response 4) mediated induction of aberrant amyloid peptide production in Alzheimer disease (AD; presenile and senile dementia). PAR4 is a leucine zipper protein that is pro-apoptotic and associated with neuronal degeneration in AD. Che-1 interaction with PAR4 suggests that it might directly or indirectly participate in regulation of PAR4 activity. Che-1 also co-localizes with PAR4 in both cytoplasmic and nuclear compartments, and it interacts directly and selectively with PAR4 via the leucine zipper domain in neural cells.

REFERENCES

- Lindfors, K., Halttunen, T., Huotari, P., Nupponen, N., Vihinen, M., Visakorpi, T., Mäki, M. and Kainulainen, H. 2000. Identification of novel transcription factor-like gene from human intestinal cells. *Biochem. Biophys. Res. Commun.* 276: 660-666.
- Di Padova, M., Bruno, T., De Nicola, F., Iezzi, S., D'Angelo, C., Gallo, R., Nicosia, D., Corbi, N., Biroccio, A., Floridi, A., Passananti, C. and Fanciulli, M. 2003. Che-1 arrests human colon carcinoma cell proliferation by displacing HDAC1 from the p21WAF1/CIP1 promoter. *J. Biol. Chem.* 278: 36496-36504.
- Xie, J. and Guo, Q. 2004. AATF protects neural cells against oxidative damage induced by amyloid β -peptide. *Neurobiol. Dis.* 16: 150-157.
- Guo, Q. and Xie, J. 2004. AATF inhibits aberrant production of amyloid β peptide 1-42 by interacting directly with PAR4. *J. Biol. Chem.* 279: 4596-4603.
- Burgdorf, S., Leister, P. and Scheidtmann, K.H. 2004. Tsg101 interacts with apoptosis-antagonizing transcription factor and enhances androgen receptor-mediated transcription by promoting its monoubiquitination. *J. Biol. Chem.* 279: 17524-17534.
- Nousiainen, M., Silljé, H.H., Sauer, G., Nigg, E.A. and Körner, R. 2006. Phosphoproteome analysis of the human mitotic spindle. *Proc. Natl. Acad. Sci. USA* 103: 5391-5396.
- Kaul, D. and Mehrotra, A. 2007. Functional characterization of AATF transcriptome in human leukemic cells. *Mol. Cell. Biochem.* 297: 215-220.

CHROMOSOMAL LOCATION

Genetic locus: AATF (human) mapping to 17q12 ; Aatf (mouse) mapping to 11 C.

SOURCE

Che-1 (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Che-1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-70032 X, 200 μ g/0.1 ml.

Blocking peptide available for competition studies, sc-70032 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Che-1 (C-19) is recommended for detection of Che-1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Che-1 (C-19) is also recommended for detection of Che-1 in additional species, including bovine, porcine and avian.

Suitable for use as control antibody for Che-1 siRNA (h): sc-72888, Che-1 siRNA (m): sc-72889, Che-1 shRNA Plasmid (h): sc-72888-SH, Che-1 shRNA Plasmid (m): sc-72889-SH, Che-1 shRNA (h) Lentiviral Particles: sc-72888-V and Che-1 shRNA (m) Lentiviral Particles: sc-72889-V.

Che-1 (C-19) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of Che-1: 63 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.


 MONOS
Satisfaction
Guaranteed

Try **Che-1 (AATF2B6): sc-81225**, our highly recommended monoclonal alternative to Che-1 (C-19).