

AKIP (C-19): sc-70106

BACKGROUND

AKIP (AURKA-interacting protein), also known as AURKAIP1 (aurora kinase A interacting protein 1) or AIP, is a 199 amino acid protein that localizes to the nucleus and is ubiquitously expressed, with highest levels present in testis, heart and skeletal muscle. Interacting specifically with ARK-1 (aurora kinase 1), AKIP functions to induce the proteasomal-dependent degradation of ARK-1, thereby acting as a negative regulator of ARK-1 activity. AKIP is encoded by a gene which maps to human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome. Chromosome 1 houses a large number of disease-associated genes, including those that are involved in familial adenomatous polyposis, Stickler syndrome, Parkinson's disease, Gaucher disease, schizophrenia and Usher syndrome. Aberrations in chromosome 1 are found in a variety of cancers, including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

- Kiat, L.S., Hui, K.M. and Gopalan, G. 2002. Aurora A kinase interacting protein (AIP), a novel negative regulator of human Aurora A kinase. *J. Biol. Chem.* 277: 45558-45565.
- Hirota, T., Kunitoku, N., Sasayama, T., Marumoto, T., Zhang, D., Nitta, M., Hatakeyama, K. and Saya, H. 2003. Aurora A and an interacting activator, the LIM protein Ajuba, are required for mitotic commitment in human cells. *Cell* 114: 585-598.
- Sun, T., Miao, X., Wang, J., Tan, W., Zhou, Y., Yu, C. and Lin, D. 2004. Functional Phe31Ile polymorphism in Aurora A and risk of breast carcinoma. *Carcinogenesis* 25: 2225-2230.
- Katayama, H., Sasai, K., Czerniak, B.A., Carter, J.L. and Sen, S. 2007. Aurora A kinase phosphorylation of Aurora A kinase interacting protein (AIP) and stabilization of the enzyme-substrate complex. *J. Cell. Biochem.* 102: 1318-1331.
- Lim, S.K. and Gopalan, G. 2007. Antizyme1 mediates AURKAIP1-dependent degradation of Aurora A. *Oncogene* 26: 6593-6603.
- Chan, E.H., Santamaria, A., Silljé, H.H. and Nigg, E.A. 2008. Plk1 regulates mitotic Aurora A function through β TrCP-dependent degradation of hBora. *Chromosoma* 117: 457-469.

CHROMOSOMAL LOCATION

Genetic locus: AURKAIP1 (human) mapping to 1p36.33; Aurkaip1 (mouse) mapping to 4 E2.

SOURCE

AKIP (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of AKIP of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-70106 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

AKIP (C-19) is recommended for detection of AKIP of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

AKIP (C-19) is also recommended for detection of AKIP in additional species, including canine and bovine.

Suitable for use as control antibody for AKIP siRNA (h): sc-72472, AKIP siRNA (m): sc-72473, AKIP shRNA Plasmid (h): sc-72472-SH, AKIP shRNA Plasmid (m): sc-72473-SH, AKIP shRNA (h) Lentiviral Particles: sc-72472-V and AKIP shRNA (m) Lentiviral Particles: sc-72473-V.


Molecular Weight of AKIP: 22 kDa.

Positive Controls: human prostate extract: sc-363774.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



AKIP (C-19): sc-70106. Western blot analysis of AKIP expression in human prostate tissue extract.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.