

Atg12 (C-18): sc-70128

BACKGROUND

Atg12 (autophagy-related protein 12), also known as APG12, APG12L, FBR93 or HAPG12, is a 140 amino acid protein that is ubiquitously expressed and belongs to the Atg12 family of proteins. Atg12 is a homolog of the yeast protein Apg12 that participates in autophagy. Autophagy is a membrane trafficking mechanism that delivers cytoplasmic cargo to the vacuole/lysosome for degradation and recycling. In yeast, autophagy requires a protein conjugation system consisting of Apg12 covalently bound at the carboxy-terminal glycine to Lysine 149 of Apg5. Similarly in humans, Atg12 is essential for autophagy and localizes to the cytoplasm where it is covalently bound to APG5, a conjugation reaction that requires APG7, Atg10 and ATP. The Atg12-APG5 conjugate functions as an important regulator of the autophagic process and is required for the change in membrane morphology and development of autophagosomes. Due to alternative splicing events, two Atg12 isoforms exist.

REFERENCES

1. Ueno, K., et al. 1998. Cloning and tissue expression of cDNAs from chromosome 5q21-22 which is frequently deleted in advanced lung cancer. *Hum. Genet.* 102: 63-68.
2. Mizushima, N., et al. 1998. A new protein conjugation system in human. The counterpart of the yeast Apg12p conjugation system essential for autophagy. *J. Biol. Chem.* 273: 33889-33892.
3. Mizushima, N., et al. 1998. A protein conjugation system essential for autophagy. *Nature* 395: 395-398.
4. Mizushima, N., et al. 2001. Dissection of autophagosome formation using Apg5-deficient mouse embryonic stem cells. *J. Cell Biol.* 152: 657-668.
5. Mizushima, N., et al. 2002. Mouse Apg10 as an Apg12-conjugating enzyme: analysis by the conjugation-mediated yeast two-hybrid method. *FEBS Lett.* 532: 450-454.
6. Tanida, I., et al. 2002. Human Apg3p/Aut1p homologue is an authentic E2 enzyme for multiple substrates, GATE-16, GABARAP, and MAP-LC3, and facilitates the conjugation of hApg12p to hApg5p. *J. Biol. Chem.* 277: 13739-13744.
7. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609608. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
8. Rabkin, S.W. and Klassen, S.S. 2007. Nitric oxide differentially regulates the gene expression of caspase genes but not some autophagic genes. *Nitric Oxide* 16: 339-347.
9. Jounai, N., et al. 2007. The Atg5 Atg12 conjugate associates with innate antiviral immune responses. *Proc. Natl. Acad. Sci. USA* 104: 14050-14055.

CHROMOSOMAL LOCATION

Genetic locus: ATG12 (human) mapping to 5q22.3; Atg12 (mouse) mapping to 18 C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

Atg12 (C-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Atg12 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-70128 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Atg12 (C-18) is recommended for detection of Atg12 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Atg12 (C-18) is also recommended for detection of Atg12 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for Atg12 siRNA (h): sc-72578, Atg12 siRNA (m): sc-72579, Atg12 shRNA Plasmid (h): sc-72578-SH, Atg12 shRNA Plasmid (m): sc-72579-SH, Atg12 shRNA (h) Lentiviral Particles: sc-72578-V and Atg12 shRNA (m) Lentiviral Particles: sc-72579-V.

Molecular Weight of Atg12 monomer: 21 kDa.

Molecular Weight of Atg12-APG5 conjugate: 60 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.



MONOS
Satisfaction
Guaranteed

Try **Atg12 (C-6): sc-271688**, our highly recommended monoclonal alternative to Atg12 (C-18).