



DACH1 (S-14): sc-70258

BACKGROUND

DACH1 (Dachshund homolog 1), also known as DACH, is a 758 amino acid homolog of the *Drosophila* dachshund gene that encodes a nuclear factor involved in eye, leg and nervous system development. Localized to the nucleus and expressed throughout the body, DACH1 is a transcription factor that regulates the activation of a variety of genes involved in organogenesis and is crucial in proper eye formation. Through association with Smad4 and NCOR1, DACH1 is able to inhibit the TGF β signaling pathway and, via its DACHbox-N domain, can bind directly to chromatin, where it regulates transcription. Additionally, DACH1 can block cellular proliferation and growth of human breast cancer cells, suggesting a possible role in tumor suppression. Four isoforms of DACH1 exist due to alternative splicing events.

REFERENCES

1. Ayres, J.A., Shum, L., Akarsu, A.N., Dashner, R., Takahashi, K., Ikura, T., Slavkin, H.C. and Nuckolls, G.H. 2001. DACH: genomic characterization, evaluation as a candidate for postaxial polydactyly type A2, and developmental expression pattern of the mouse homologue. *Genomics* 77: 18-26.
2. Heanue, T.A., Davis, R.J., Rowitch, D.H., Kispert, A., McMahon, A.P., Mardon, G. and Tabin, C.J. 2002. DACH1, a vertebrate homologue of *Drosophila* dachshund, is expressed in the developing eye and ear of both chick and mouse and is regulated independently of Pax and Eya genes. *Mech. Dev.* 111: 75-87.
3. Wu, K., Yang, Y., Wang, C., Davoli, M.A., D'Amico, M., Li, A., Cveklova, K., Kozmik, Z., Lisanti, M.P., Russell, R.G., Cvekl, A. and Pestell, R.G. 2003. DACH1 inhibits transforming growth factor β signaling through binding Smad4. *J. Biol. Chem.* 278: 51673-51684.
4. Wu, K., Li, A., Rao, M., Liu, M., Dailey, V., Yang, Y., Di Vizio, D., Wang, C., Lisanti, M.P., Sauter, G., Russell, R.G., Cvekl, A. and Pestell, R.G. 2006. DACH1 is a cell fate determination factor that inhibits cyclin D1 and breast tumor growth. *Mol. Cell. Biol.* 26: 7116-7129.
5. Sunde, J.S., Donninger, H., Wu, K., Johnson, M.E., Pestell, R.G., Rose, G.S., Mok, S.C., Brady, J., Bonome, T. and Birrer, M.J. 2006. Expression profiling identifies altered expression of genes that contribute to the inhibition of transforming growth factor β signaling in ovarian cancer. *Cancer Res.* 66: 8404-8412.

CHROMOSOMAL LOCATION

Genetic locus: DACH1 (human) mapping to 13q22; Dach1 (mouse) mapping to 14 E3.

SOURCE

DACH1 (S-14) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of DACH1 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-70258 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-70258 X, 200 μ g/0.1 ml.

APPLICATIONS

DACH1 (S-14) is recommended for detection of DACH1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DACH1 siRNA (h): sc-77089, DACH1 siRNA (m): sc-77090, DACH1 shRNA Plasmid (h): sc-77089-SH, DACH1 shRNA Plasmid (m): sc-77090-SH, DACH1 shRNA (h) Lentiviral Particles: sc-77089-V and DACH1 shRNA (m) Lentiviral Particles: sc-77090-V.

DACH1 (S-14) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of DACH1: 79 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.