

GPR4 (S-14): sc-70290

BACKGROUND

G protein-coupled receptors (GPRs), also known as seven transmembrane receptors, heptahelical receptors or 7TM receptors, comprise a superfamily of proteins that play a role in many different stimulus-response pathways. G protein coupled receptors translate extracellular signals into intracellular signals (G protein activation) and they respond to a variety of signaling molecules, such as hormones and neurotransmitters. GPR4 (G protein-coupled receptor 4), also known as GPR19, is a 362 amino acid multi-pass membrane protein that belongs to the G protein-coupled receptor family and functions as an orphan receptor. GPR4 is overexpressed in cancer tissue, suggesting a role in tumor formation, and is also believed to play a critical role in endothelial cell function and may mediate the effects of sphingosylphosphorylcholine.

REFERENCES

1. An, S., Tsai, C. and Goetzl, E.J. 1995. Cloning, sequencing and tissue distribution of two related G protein-coupled receptor candidates expressed prominently in human lung tissue. *FEBS Lett.* 375: 121-124.
2. Mahadevan, M.S., Baird, S., Bailly, J.E., Shutler, G.G., Sabourin, L.A., Tsilfidis, C., Neville, C.E., Narang, M. and Korneluk, R.G. 1995. Isolation of a novel G protein-coupled receptor (GPR4) localized to chromosome 19q13.3. *Genomics* 30: 84-88.
3. Zhu, K., Baudhuin, L.M., Hong, G., Williams, F.S., Cristina, K.L., Kabarowski, J.H., Witte, O.N. and Xu, Y. 2001. Sphingosylphosphorylcholine and lysophosphatidylcholine are ligands for the G protein-coupled receptor GPR4. *J. Biol. Chem.* 276: 41325-41335.
4. Bektas, M., Barak, L.S., Jolly, P.S., Liu, H., Lynch, K.R., Lacana, E., Suhr, K.B., Milstien, S. and Spiegel, S. 2003. The G protein-coupled receptor GPR4 suppresses ERK activation in a ligand-independent manner. *Biochemistry* 42: 12181-12191.
5. Online Mendelian Inheritance in Man, OMIM™. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 600551. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Sin, W.C., Zhang, Y., Zhong, W., Adhikarakunnathu, S., Powers, S., Hoey, T., An, S. and Yang, J. 2004. G protein-coupled receptors GPR4 and TDAG8 are oncogenic and overexpressed in human cancers. *Oncogene* 23: 6299-6303.
7. Kim, K.S., Ren, J., Jiang, Y., Ebrahem, Q., Tipps, R., Cristina, K., Xiao, Y.J., Qiao, J., Taylor, K.L., Lum, H., Anand-Apte, B. and Xu, Y. 2005. GPR4 plays a critical role in endothelial cell function and mediates the effects of sphingosylphosphorylcholine. *FASEB J.* 19: 819-821.
8. Qiao, J., Huang, F., Naikawadi, R.P., Kim, K.S., Said, T. and Lum, H. 2006. Lysophosphatidylcholine impairs endothelial barrier function through the G protein-coupled receptor GPR4. *Am. J. Physiol. Lung Cell. Mol. Physiol.* 291: L91-L101.

CHROMOSOMAL LOCATION

Genetic locus: GPR4 (human) mapping to 19q13.32; Gpr4 (mouse) mapping to 7 A3.

SOURCE

GPR4 (S-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of GPR4 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-70290 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GPR4 (S-14) is recommended for detection of GPR4 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

GPR4 (S-14) is also recommended for detection of GPR4 in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for GPR4 siRNA (h): sc-75179, GPR4 siRNA (m): sc-75180, GPR4 shRNA Plasmid (h): sc-75179-SH, GPR4 shRNA Plasmid (m): sc-75180-SH, GPR4 shRNA (h) Lentiviral Particles: sc-75179-V and GPR4 shRNA (m) Lentiviral Particles: sc-75180-V.

Molecular Weight of GPR4: 45 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.