



Blood Group M antigen (1.B.710): sc-70428

BACKGROUND

Blood-group antigens are generally defined as molecules formed by sequential addition of saccharides to the carbohydrate side chains of lipids and proteins detected on erythrocytes and certain epithelial cells. The A, B and H antigens are reported to undergo modulation during malignant cellular transformation. Blood group related antigens are usually mucin-type and are detected on erythrocytes, certain epithelial cells and in secretions of certain individuals. Sixteen genetically and biosynthetically distinct but inter-related specificities belong to this group of antigens, including A (1 and 2), B, H, M, N, Lewis A, Lewis B, Lewis X, Lewis Y and precursor type 1 chain antigens.

REFERENCES

1. Watkins, W.M. 1980. Biochemistry and genetics of the ABO, Lewis: P blood group systems. *Adv. Hum. Genet.* 10: 379-385.
2. Fraser, R.H., Inglis, G., Mackie, A., Munro, A.C., Allan, E.K., Mitchell, R., Sonneborn, H.H. and Uthemann, H. 1985. Mouse monoclonal antibodies reacting with M blood group-related antigens. *Transfusion* 25: 261-266.
3. Jaskiewicz, E., Moulds, J.J., Kraemer, K., Goldstein, A.S. and Lisowska, E. 1990. Characterization of the epitope recognized by a monoclonal antibody highly specific for Blood Group M antigen. *Transfusion* 30: 230-235.
4. Wimmer, B. 1992. Controlled percutaneous puncture of gastrointestinal and hepatopancreatic lesions. *Bildgebung* 59: 68-71.
5. Boren, T., Falk, P., Roth, K.A., Larson, G. and Normark, S. 1994. Attachment of *Helicobacter pylori* to human group antigens. *Science* 262: 1892-1895.
6. Appelmek, B.J., Simoons-Smit, I., Negrini, R., Moran, A.P., Aspinall, G.O., Forte, J.G., De Vries, T., Quan, H., Verboom, T., Maaskant, J.J., Ghiara, P., Kuipers, E.J., Bloemena, E., Tadema, T.M., Townsend, R.R., Tyagarajan, K., Crothers, J.M., Monteiro, M.A., Savio, A. and De Graaff, J. 1996. Potential role of molecular mimicry between *Helicobacter pylori* lipopolysaccharide and host Lewis blood group antigens in autoimmunity. *Infect. Immun.* 64: 2031-2040.
7. Ilver, D., Arnqvist, A., Ogren, J., Frick, I.M., Kersulyte, D., Incecik, E.T., Berg, D.E., Covacci, A., Engstrand, L. and Boren, T. 1998. *Helicobacter pylori* adhesin binding fucosylated histo-blood group antigens revealed by retagging. *Science* 279: 373-377.

CHROMOSOMAL LOCATION

Genetic locus: GYPA/GYPB (human) mapping to 4q31.21.

SOURCE

Blood Group M antigen (1.B.710) is a mouse monoclonal antibody raised against M antigen on cells of human origin.

PRODUCT

Each vial contains 100 µg IgM in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Blood Group M antigen (1.B.710) is recommended for detection of Blood Group M antigen of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and flow cytometry (1 µg per 1 x 10⁶ cells).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.