

Cdk6 (H-230): sc-7181

BACKGROUND

Cell cycle progression is controlled in part by a family of cyclin proteins and cyclin dependent kinases (Cdk). Cdk proteins work in concert with the cyclins to phosphorylate key substrates involved in each phase of cell cycle progression. Another family of proteins, Cdk inhibitors, also plays a role in regulating the cell cycle by binding to cyclin/Cdk complexes and modulating their activity. Several Cdk proteins have been identified, including Cdk2-Cdk8, PCTAIRE-1-3, PITSLRE and PITSLRE. Cdk6 is known to associate with cyclins D1, D2 and D3 and to be involved with the G₁/S transition of the cell cycle. Multiple inhibitors of Cdk6 have been identified, including p18 and p19. These inhibitors bind to both free and complexed Cdk6, and they inhibit the activity of the cyclin D-bound Cdk6.

CHROMOSOMAL LOCATION

Genetic locus: CDK6 (human) mapping to 7q21.2; Cdk6 (mouse) mapping to 5 A1.

SOURCE

Cdk6 (H-230) is a rabbit polyclonal antibody raised against amino acids 1-230 mapping at the N-terminus of Cdk6 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Cdk6 (H-230) is recommended for detection of Cdk6 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Cdk6 (H-230) is also recommended for detection of Cdk6 in additional species, including equine, canine, bovine and avian.

Suitable for use as control antibody for Cdk6 siRNA (h): sc-29264, Cdk6 siRNA (m): sc-35048, Cdk6 shRNA Plasmid (h): sc-29264-SH, Cdk6 shRNA Plasmid (m): sc-35048-SH, Cdk6 shRNA (h) Lentiviral Particles: sc-29264-V and Cdk6 shRNA (m) Lentiviral Particles: sc-35048-V.

Molecular Weight of Cdk6: 40 kDa.

Positive Controls: K-562 whole cell lysate: sc-2203, Jurkat whole cell lysate: sc-2204 or Jurkat + PMA nuclear extract: sc-2133.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SELECT PRODUCT CITATIONS

1. Belyavskiy, M., et al. 1998. The structural protein ODV-EC27 of *Autographa californica* nucleopolyhedrovirus is a multifunctional viral cyclin. *Proc. Natl. Acad. Sci. USA* 95: 11205-11210.
2. Crish, J.F., et al. 2000. Suprabasal expression of the human papillomavirus type 16 oncoproteins in mouse epidermis alters expression of cell cycle regulatory proteins. *Carcinogenesis* 21: 1031-1037.
3. Masaki, T., et al. 2000. Hepatocellular carcinoma cell cycle: study of Long-Evans cinnamon rats. *Hepatology* 32: 711-720.
4. Bhunia, A.K., et al. 2002. PKD1 induces p21^{waf1} and regulation of the cell cycle via direct activation of the JAK/Stat signaling pathway in a process requiring PKD2. *Cell* 109: 157-168.
5. Iida, T., et al. 2002. Hypoxia-inducible factor-1 α induces cell cycle arrest of endothelial cells. *Genes Cells* 7: 143-149.
6. Nguyen, M.D., et al. 2003. Cell cycle regulators in the neuronal death pathway of amyotrophic lateral sclerosis caused by mutant superoxide dismutase 1. *J. Neurosci.* 23: 2131-2340.
7. Sah, J.F., et al. 2004. Epigallocatechin-3-gallate inhibits epidermal growth factor receptor signaling pathway: evidence for direct inhibition of ERK 1/2 and Akt kinases. *J. Biol. Chem.* 279: 12755-12762.
8. Blaschke, F., et al. 2004. Liver X receptor agonists suppress vascular smooth muscle cell proliferation and inhibit neointima formation in balloon-injured rat carotid arteries. *Circ. Res.* 95: e110-e123.
9. Balasubramanian, S., et al. 2004. Suppression of human pancreatic cancer cell proliferation by AGN194204, an RXR-selective retinoid. *Carcinogenesis* 25: 1377-1385.
10. Balasubramanian, S., et al. 2005. A novel retinoid-related molecule inhibits pancreatic cancer cell proliferation by a retinoid receptor independent mechanism via suppression of cell cycle regulatory protein function and induction of caspase-associated apoptosis. *Oncogene* 24: 4257-4270.
11. Wu, W., et al. 2009. Antibody array analysis with label-based detection and resolution of protein size. *Mol. Cell. Proteomics* 8: 245-257.
12. Baranda-Avila, N., et al. 2009. Molecular mechanism of cell proliferation in rodent uterus during the estrous cycle. *J. Steroid Biochem. Mol. Biol.* 113: 259-268.
13. Ikner, A., et al. 2011. TWEAK induces apoptosis through a death-signaling complex comprising receptor-interacting protein 1 (RIP1), Fas-associated death domain (FADD), and caspase-8. *J. Biol. Chem.* 286: 21546-21554.

MONOS
Satisfaction
Guaranteed

Try **Cdk6 (B-10): sc-7961** or **Cdk6 (DCS-83): sc-53638**, our highly recommended monoclonal alternatives to Cdk6 (H-230). Also, for AC, HRP, FITC, PE, Alexa Fluor[®] 488 and Alexa Fluor[®] 647 conjugates, see **Cdk6 (B-10): sc-7961**.