

p53 (3H2820): sc-71821

BACKGROUND

p53, a DNA-binding, oligomerization domain- and transcription activation domain-containing tumor suppressor that upregulates growth arrest and apoptosis-related genes in response to stress signals, thereby influencing programmed cell death, cell differentiation and cell cycle control mechanisms. p53 localizes to the nucleus, yet can be chaperoned to the cytoplasm by the negative regulator MDM2, an E3 ubiquitin ligase that is upregulated in the presence of active p53, where MDM2 polyubiquitinates p53 for proteasome targeting. p53 fluctuates between latent and active (DNA-binding) conformations, and is differentially activated through post-translational modifications including phosphorylation and acetylation. Mutations in the DNA-binding domain (DBD) of p53, amino acids 110-286, can compromise energetically favorable association with *cis* elements and are implicated in several human cancers.

REFERENCES

1. Banks, L., et al. 1986. Isolation of human p53-specific monoclonal antibodies and their use in the studies of human p53 expression. *Eur. J. Biochem.* 159: 529-534.
2. Hupp, T.R., et al. 1992. Regulation of the specific DNA binding function of p53. *Cell* 71: 875-886.
3. Levine, A.J. 1997. P53, the cellular gatekeeper for growth and division. *Cell* 88: 323-331.
4. Ashcroft, M., et al. 1999. Regulation of p53 stability. *Oncogene* 18: 7637-7643.
5. Soussi, T., et al. 2000. P53 website and analysis of p53 gene mutations in human cancer: forging a link between epidemiology and carcinogenesis. *Hum. Mutat.* 15: 105-113.
6. Chene, P. 2001. The role of tetramerization in p53 function. *Oncogene* 20: 2611-2617.
7. Minamoto, T., et al. 2001. Distinct pattern of p53 phosphorylation in human tumors. *Oncogene* 20: 3341-3347.
8. LocusLink Report (LocusID: 7157). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: TP53 (human) mapping to 17p13.1; Trp53 (mouse) mapping to 11 B3.

SOURCE

p53 (3H2820) is a mouse monoclonal antibody raised against p53- β galactosidase fusion protein.

PRODUCT

Each vial contains 200 μ g IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS


p53 (3H2820) is recommended for detection of p53 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500).

Suitable for use as control antibody for p53 siRNA (h): sc-29435, p53 siRNA (m): sc-29436, p53 siRNA (r): sc-45917, p53 shRNA Plasmid (h): sc-29435-SH, p53 shRNA Plasmid (m): sc-29436-SH, p53 shRNA Plasmid (r): sc-45917-SH, p53 shRNA (h) Lentiviral Particles: sc-29435-V, p53 shRNA (m) Lentiviral Particles: sc-29436-V and p53 shRNA (r) Lentiviral Particles: sc-45917-V.


Molecular Weight of p53: 53 kDa.

Positive Controls: mouse LacZ whole cell lysate: sc-364371, A-431 whole cell lysate: sc-2201 or p53 (m): 293T Lysate: sc-125766.

DATA



p53 (3H2820): sc-71821. Western blot analysis of p53 expression in A-431 (A) and mouse LacZ (B) whole cell lysates.



p53 (3H2820): sc-71821. Western blot analysis of p53 expression in non-transfected: sc-117752 (A) and mouse p53 transfected: sc-125766 (B) 293T whole cell lysates.

SELECT PRODUCT CITATIONS

1. Dhar, S., et al. 2016. Dietary pterostilbene is a novel MTA1-targeted chemopreventive and therapeutic agent in prostate cancer. *Oncotarget* 7: 18469-18484.
2. Bastos, C.C.C., et al. 2016. Use of *Bidens pilosa* L. (Asteraceae) and *Curcuma longa* L. (Zingiberaceae) to treat intestinal mucositis in mice: toxicopharmacological evaluations. *Toxicol. Rep.* 3: 279-287.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CONJUGATES

See **p53 (DO-1): sc-126** for p53 antibody conjugates, including AC, HRP, FITC, PE, and Alexa Fluor® 488, 546, 594, 647, 680 and 790.