

Factor I siRNA (m): sc-72878

BACKGROUND

The complement pathway is an important host defense system that contributes to both innate and acquired immunity. There are three pathways of complement activation: the classical pathway, lectin pathway and alternative pathway. Complement protein Factor I is a key serine protease that modulates the complement cascade by regulating the levels of C3 convertases. It circulates in plasma as a heavily N-glycosylated heterodimer made up of two disulfide linked chains, each carrying three N-linked oligosaccharide chains that may have both structural and functional roles in the interactions with the natural substrate and the cofactor during catalysis. Factor I is a serine protease with a high degree of specificity for C3b and C4b. It requires protein cofactors for cleavage of these complement proteins; Factor H, CR1 or MCP are required for C3b cleavage, and C4bp or CR1 are required for C4b cleavage.

REFERENCES

- Schlaf, G., et al. 2001. Complement Factor I is upregulated in rat hepatocytes by interleukin-6 but not by interferon- γ , interleukin-1 β , or tumor necrosis factor- α . *Biol. Chem.* 382: 1089-1094.
- Terado, T., et al. 2002. Conservation of the modular structure of complement Factor I through vertebrate evolution. *Dev. Comp. Immunol.* 26: 403-413.
- Cunnion, K.M., et al. 2004. Cleavage of complement C3b to iC3b on the surface of *Staphylococcus aureus* is mediated by serum complement Factor I. *Infect. Immun.* 72: 2858-2863.
- Fremeaux-Bacchi, V., et al. 2004. Complement Factor I: a susceptibility gene for atypical haemolytic uraemic syndrome. *J. Med. Genet.* 41: e84.
- Genel, F., et al. 2005. Complement Factor I and immune complex glomerulonephritis. *Scand. J. Infect. Dis.* 37: 615-618.
- Tsiftoglou, S.A., et al. 2005. The catalytically active serine protease domain of human complement Factor I. *Biochemistry* 44: 6239-6249.
- Grumach, A.S., et al. 2006. Recurrent infections in partial complement factor I deficiency: evaluation of three generations of a Brazilian family. *Clin. Exp. Immunol.* 143: 297-304.
- Liu, W., et al. 2006. Down-regulated expression of complement factor I: a potential suppressive protein for gastric cancer identified by serum proteome analysis. *Clin. Chim. Acta* 377: 119-126.
- Tsiftoglou, S.A., et al. 2006. Human complement factor I glycosylation: structural and functional characterisation of the N-linked oligosaccharides. *Biochim. Biophys. Acta* 1764: 1757-1766.

CHROMOSOMAL LOCATION

Genetic locus: Cfi (mouse) mapping to 3 G3.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Factor I siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Factor I shRNA Plasmid (m): sc-72878-SH and Factor I shRNA (m) Lentiviral Particles: sc-72878-V as alternate gene silencing products.

For independent verification of Factor I (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-72878A, sc-72878B and sc-72878C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Factor I siRNA (m) is recommended for the inhibition of Factor I expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Factor I gene expression knockdown using RT-PCR Primer: Factor I (m)-PR: sc-72878-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.