



Laminin-5 (B5815M): sc-73626

BACKGROUND

Laminin-5 is a glycoprotein complex of three subunits (Laminin α -3, β -3, and γ -2) that influences cell adhesion (metastasis), signal transduction and keratinocyte differentiation. Laminin-5 localizes to the basal lamina underneath epithelia and mediates the anchoring of basal epithelial cells to the extracellular matrix (ECM). Differential processing of the subunits of the Laminin-5 precursor influences how this protein integrates into the ECM architecture.

REFERENCES

- Giannelli, G., Fransvea, E., Bergamini, C., Marinosci, F. and Antonaci, S. 2003. Laminin-5 chains are expressed differentially in metastatic and non-metastatic hepatocellular carcinoma. *Clin. Cancer Res.* 9: 3684-3691.
- Tsuruta, D., Hopkinson, S.B., Lane, K.D., Werner, M.E., Cryns, V.L. and Jones, J.C. 2003. Crucial role of the specificity-determining loop of the Integrin β 4 subunit in the binding of cells to Laminin-5 and outside-in signal transduction. *J. Biol. Chem.* 278: 38707-38714.
- Aumailley, M., El Khal, A., Knoss, N. and Tunggal, L. 2003. Laminin-5 processing and its integration into the ECM. *Matrix Biol.* 22: 49-54.
- Elkhal, A., Tunggal, L. and Aumailley, M. 2004. Fibroblasts contribute to the deposition of Laminin-5 in the extracellular matrix. *Exp. Cell Res.* 296: 223-230.
- Kariya, Y. and Miyazaki, K. 2004. The basement membrane protein Laminin-5 acts as a soluble cell motility factor. *Exp. Cell Res.* 297: 508-520.
- Wang, H., Fu, W., Im, J.H., Zhou, Z., Santoro, S.A., Iyer, V., DiPersio, C.M., Yu, Q.C., Quaranta, V., Al-Mehdi, A. and Muschel, R.J. 2004. Tumor cell α 3/ β 1 Integrin and vascular Laminin-5 mediate pulmonary arrest and metastasis. *J. Cell Biol.* 164: 935-941.
- Kunneken, K., Pohlentz, G., Schmidt-Hederich, A., Odenthal, U., Smyth, N., Peter-Katalinic, J., Bruckner, P. and Eble, J.A. 2004. Recombinant human Laminin-5 domains. Effects of heterotrimerization, proteolytic processing, and N-glycosylation on α 3/ β 1 Integrin binding. *J. Biol. Chem.* 279: 5184-5193.

CHROMOSOMAL LOCATION

Genetic locus: LAMA5 (human) mapping to 20q13.2-q13.3.

SOURCE

Laminin-5 (B5815M) is a mouse monoclonal antibody raised against Laminin-5 β 3 chain of human origin.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

Laminin-5 (B5815M) is recommended for detection of Laminin-5 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000).

Suitable for use as control antibody for Laminin-5 siRNA (h): sc-35788, Laminin-5 shRNA Plasmid (h): sc-35788-SH and Laminin-5 shRNA (h) Lentiviral Particles: sc-35788-V.

Molecular Weight of Laminin-5: 170 kDa.

Positive Controls: A-431 whole cell lysate: sc-2201.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.